LED Digital 7" Square Hotplate Stirrer

User Manual



MS7-S	Classic 7" Square Plate Stirrer
MS7-H550-S	LED Digital 7" Square Hotplate Stirrer

Please read the User Manual carefully before use, and follow all operating and safety instructions!

www.scilogex.com



Contents

Pref	ace1
Serv	ice1
War	ranty1
1.	Safety Instructions
2.	Proper use
3.	Inspection4
	3.1 Receiving Inspection4
	3.2 Listing of Items
4.	Control
	4.1 Control elements
	4.2 Display (Hotplate Model)6
5.	Trial run6
6.	Working with external temperature sensor (Hotplate
Moc	lel)7
7.	Faults
8.	Maintenance and Cleaning8
9.	Associated standards and regulations9
10.	Technical data9
11.	Ordering information10

LED Digital 7" Square Hotplate Stirrer

Preface

Welcome to the "LED Digital 7" Square Hotplate Stirrer". Users should read this Manual carefully, follow the instructions and procedures, and be aware of all the cautions when using this instrument.

Service

When help needed, you can always contact the Service Department of manufacturer for technical support in the following ways:

SCILOGEX, LLC 500 Four Rod Road Suite 122 Berlin, CT 06037 USA Tel: 1- (860) 828-5614 Fax: 1- (860) 828-5389 E-mail: info@scilogex.com Website: Http://www.scilogex.com Please provide the customer care representative with the following information:

- Serial number (on the rear panel)
- Certification
- Description of problem (i.e., hardware or software)
- Methods and procedures adopted to resolve the problems
- Your contact information

Warranty

This instrument is warranted to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation. For claims under the warranty please contact your local supplier. You may also send the instrument directly to manufacturer, enclosing the invoice copy and by giving reasons for the claim.

1. Safety Instructions

Warning!

- Read the operating instructions carefully before use.
- Ensure that only trained staff works with the instrument.

Risk of burn!



- Caution when touch the housing parts and the hotplate which can reach temperature of 550 ℃.
- Pay attention to the residual heat after switching off.

Protective ground contact!



• Make sure that socket must be grounded (protective ground contact) before use.

• When working wear personal safety guards to avoid the risk from:

- Splashing and evaporation of liquids
- Release of toxic or combustible gases

• Set up the instrument in a spacious are on a stable, clean, non-slip, dry and fireproof surface. Do not operate the instrument in explosive atmospheres, with hazardous substances or under water.

- Gradually increase the speed, reduce the speed if:
 - Stirring bar breaks away due to high speed
- The instrument is not running smoothly, or container moves on the base plate
- Temperature must always be set to at least 50° C lower than the fire point of the media used.
- Be aware of hazards due to:
- Flammable materials or media with a low boiling temperature
 - Overfilling of media
- Unsafe container
- Process pathogenic materials only in closed vessels.
- Check the instrument and accessories prior to each use. Do not use damaged components. Safe operation is only guaranteed with the accessories described in the "Accessories" chapter. Accessories must be securely attached to the device and can not come off by themselves. Always disconnect the plug before fitting accessories.
- When the external temperature sensor needed, the tip of the measuring sensor must be at least 5-10mm from vessel bottom and wall.

- The instrument can only be disconnected from the main power supply by pulling out the main or the connector plug.
- The voltage stated on the label must correspond to the main power supply.
- Ensure that the main power supply cable does not touch the hotplate. Do not cover the device.
- Forbid to put pressure and over heat media on the surface of glass ceramic, that can be caused surface broken.
- The instrument may only be opened by experts.
- Keep away from high magnetic field.

2. Proper use

The instrument is designed for mixing and/or heating liquids in schools, laboratories or factories.

• Observe the minimum distances between the devices, between the device and the wall and above the assembly (min. 100 mm)




Figure 1

This device is not suitable for using in residential areas or other constraints mentioned in Chapter 1.

3. Inspection

3.1 Receiving Inspection

Unpack the equipment carefully and check for any damages which may have arisen during transport. Please contact manufacturer/supplier for technical support.



Note:

If there is any apparent damage to the system, please do not plug it into the power line.

3.2 Listing of Items

The package includes the following items:

Items	Qty	
Main unit	1	
Power cable	1	
User Manual	1	
Table 1		


4. Control

4.1 Control elements



Hotplate model

Figure 2






Figure 5

Items	Descriptions	
Speed control Knob	The stirring function is switched	
Stir	ON or OFF by rotating the knob.	
Temperature control	The heating function is switched	
knob Heat (Hotplate)	ON or OFF by rotating the knob.	
	If rotate the heating knob, LED	
	displays the temperature setting	
LED (Hotplate)	value and shift to real value in the	
	duration of 5 seconds.	
	When the heating function is	
LED Heat (Hotplate)	switched ON, the LED Heat	
	flash.	

	When the instrument is
LED Stir	switched ON, the LED Stir is
	lit.
Deserve Carritel	Switch ON or OFF the
Power Switch	instrument.
	Table 2

Table 2

4.2 Display (Hotplate Model)



Figure 6

Display	Descriptions		
Display area	When the heating function is switched ON and rotate the stirring knob, LED displays the temperature setting value and shifts to real value in 5 seconds. When the heating function is switched OFF and the hotplate temperature is still above 50°C, LED displays Hot, otherwise LED displays 0.		

5. Trial run

• Make sure the required operating voltage and power supply voltage match.

• Ensure the socket must be properly grounded.

• Plug in the power cable, ensure the power is on and begin initializing.

• Add the medium into the vessel with an appropriate stirring bar.

• Place vessel on the work plate.

• Clockwise rotate the stirring knob slowly to set the rated rotary speed in the safe speed range from 0 to 1500rpm, and start stirring function.

• Clockwise rotate the heating knob slowly to set the rated temperature in the safe range from room temperature to 550° C, and start heating function (hotplate model).

• Begin the heating and stirring functions. If these operations above are normal, the device is ready to operate. If these operations are not normal, the device may be damaged during transportation, please contact manufacturer/supplier for technical support.



Warning !

Forbid to transfer the vessel when the instrument working.

6. Working with external temperature sensor (Hotplate Model)




Figure 7

The external temperature sensor PT1000 is the Manufacture's standard accessory. If the sensor is plugged in and rotate the heating knob, LED displays the temperature setting value and shifts to real value in 5 seconds. Safe circuit controls hotplate temperature. Comparing with the temperature control of the hotplate, the external temperature sensor can control the medium's temperature more precisely.

If external temperature sensors needed, it must be plugged in before the device is switched ON. Ensure the external temperature sensors connecting firmly and to be inserted in the media before heating. The heating function will be stopped automatically under abnormal conditions. Please operate follow the instructions below:

- Switch OFF the instrument.
- Ensure the external temperature sensors is inserted in the media heated.

• Switch ON the instrument and run heating function. If the heating function did not work, please contact manufacturer/supplier for technical support.

Warning !



Forbid to plug in/pull out external temperature sensor in heating.

7. Faults

- Instruments can not be power ON.
 - Check whether the power cable is plugged in.
 - Check whether the fuse is broken or loose.

If these faults are not resolved, please contact manufacturer/supplier.

8. Maintenance and Cleaning

- Proper maintenance can keep instruments working properly and lengthen its lifetime.
- Do not spray cleanser into the instrument when cleaning.
- Unplug the power line when cleaning.
- Only use recommended cleansers:

Dyes	Isopropyl alcohol	
Construction materials	Water containing tenside	
	/ Isopropyl alcohol	
Cosmetics	Water containing tenside	
	/ Isopropyl alcohol	
Foodstuffs	Water containing tenside	
Fuels	Water containing tenside	

• Wear the proper protective gloves during cleaning of

the instrument.

- Before using other method for cleaning or decontamination, the user must ascertain with the manufacturer that this method will not harm the instrument.
- Send in the case of service the instrument back in the packaging carton. Storage packing is not sufficient for the back dispatch. Use additionally a suitable transportation packing.
- The enamel makes the hotplate easier to care for and more resistant to acids and bases. Because of it, however, the hotplate is also more susceptible to extreme fluctuations in temperature and the force of impact. This can result in cracks forming or the coating flaking off.

Warning!

Cut off power when maintenance and cleaning.

9. Associated standards and	10.Technical data		
regulations	Voltage [VAC]	100-120/200-240	
Construction in accordance with the following safety	Frequency [Hz]	50/60	
standards:	Deserver (W/1	*1030	
EN 61010-1	Power [W]	30	
UL 3101-1	Stirring point position quantity	1	
CAN/CSA C22.2(1010-1)	Max. stirring quantity (H ₂ O) [l]	10	
EN 61010-2-10 Construction in accordance with the following EMC	Max. magnetic bar [Lר, mm]	80×10	
standards:	Motor type	Shaded pole Motor	
EN 61326-1	Max. power input of motor [W]	15	
Associated EU guidelines:	Max. power output of motor [W]	1.5	
EMC-guidelines: 89/336/EWG	Speed range [rpm]	0-1500	
Instrument guidelines: 73/023/EWG	Rotary speed display	Scale	
	Plate material	Glass ceramic	
	Dimensions of workplate (mm)	184×184	

*Heating power [W]

*Temperature range [\mathbb{C}]

*Temperature display [$^{\circ}$ C]

*Temperature display accuracy [$^{\circ}$ C]

1000

LED

 ± 1

RT-550

*The safety temperature of the	580 11. Order		ing information	
hotplate [°C]	580			
*Temperature sensor in medium	PT1000	Cat No.	Descriptions	
*Control accuracy of heating			MS7-H550-S LED Digital 7'' Square	
temperature with temperature	± 0.5	813221009999	Hotplate Magnetic Stirrer, glass ceramics	
sensor [°C]			hotplate, USA plug, 110V/50Hz/60Hz MS7-H550-S LED Digital 7'' Square	
*Residual heat warning	50°C	813221019999	Hotplate Magnetic Stirrer, glass ceramics	
Dimensions (mm)	215×360×112	013221017777	hotplate, Cn plug, 220V/50Hz/60Hz	
	*4.5		MS7-H550-S LED Digital 7" Square	
Weight [kg]	3.8 813221029999	813221029999	Hotplate Magnetic Stirrer, glass ceramics	
Permitted ambient temperature [$^{\circ}$ C]	5-40		hotplate, Euro plug, 220V/50Hz/60Hz	
Permitted relative humidity	80%		MS7-H550-S LED Digital 7" Square	
Protection class acc. to DIN 60529	IP21	813221039999	Hotplate Magnetic Stirrer, glass ceramics	
		<u> </u>	hotplate, UK plug, 220V/50Hz/60Hz	
•	*Hotplate model		MS7-S Classic 7" Square Plate	
Table 3		813211009999	Magnetic Stirrer, glass ceramics plate,	
			USA plug, 110V/50Hz/60Hz	
			MS7-S Classic 7" Square Plate	
		813211019999	Magnetic Stirrer, glass ceramics plate,	
			Cn plug, 220V/50Hz/60Hz	
		813211029999	MS7-S Classic 7" Square Plate	
		013211023333	Magnetic Stirrer, glass ceramics plate,	

LED Digital 7" Square Hotplate Stirrer

Cat No.	Descriptions	Cat No.	Descriptions
	Euro plug, 220V/50Hz/60Hz		reaction vessel, Ø17.75mm, 26mm
	MS7-S Classic 7" Square Plate		depth
813211039999	Magnetic Stirrer, glass ceramics plate,		MS135.7 Golden quarter pie, 4 holes,
	UK plug, 220V/50Hz/60Hz	18900049	16ml reaction vessel, Ø21.6mm, 31.7mm
Accessories			depth
	PT1000 Temperature sensor for digital		Table 4
18900016	hotplate model, used for Pro,		
	MS7-H550-Pro & MS7-H550-S, 230mm		
	PT1000 Temperature sensor for digital		
18900084	hotplate model, used for Pro,		
	MS7-H550-Pro & MS7-H550-S, 100mm		
18900017	Support clamp of PT1000		
18900002MS135.2 Red quarter pie, 11 holes, 4 ml reaction vessel, Ø15.2mm, 20mm depth			
18000003	MS135.3 Purple quarter pie, 4 holes, 20		
18900003 ml reaction vessel, Ø28mm, 24mm depth			
MS135.4 Blue quarter pie, 4 holes, 30 ml			
18900004	18900004 reaction vessel, Ø28mm, 30mm depth		
18900005	MS135.5 Black quarter pie, 4 holes,40		
18900005	ml reaction vessel, Ø28mm, 43mm depth		
18900048	MS135.6 Green quarter pie, 6 holes, 8ml		

SCILOGEX, LLC

500 Four Rod Road Suite 122 Berlin, CT 06037 USA Tel: 1-860-828-5289 Fax:1-860-828-5389 E-mail: info@dragonlabusa.com Website: Http://www.scilogex.com