

CSF
RANGE




Contact Shock Freezers

B Medical Systems | CSF Range

Contact shock freezers are devices intended for the quick-freezing and subsequent temporary storage of blood plasma or biological samples to a core temperature below -30°C within less than one hour.

Compliant to AABB Standards | Medical Device according to 21CFR Part 864.9700, Class II



medical
systems

SAVING LIVES
THROUGH RELIABLE
AND INNOVATIVE
TECHNOLOGY

Blood Management Solutions





Safety Standards B Medical Systems



The Safety Standards developed by B Medical Systems define certain significant technical features of a product. These ensure the safe storage of the preparations as well as setting the highest standards of safety for the user.

CSF RANGE	
<input checked="" type="checkbox"/>	B Medical Systems Electronics
<input checked="" type="checkbox"/>	B Medical Systems 7" full touchscreen display
<input checked="" type="checkbox"/>	Power indicator light and digital temperature indicator (display: 0.1 digits)
<input checked="" type="checkbox"/>	Self-contained alarm system with integrated battery takes over the alarm function and temperature value measurements in case of system failure for at least 48 hours
<input checked="" type="checkbox"/>	Acoustic and visual alarm signal in case of temperature alarm and system failure
<input checked="" type="checkbox"/>	The alarm history function on the electronic stores all the relevant values during a temperature alarm, such as: min., max. and average temperature
<input checked="" type="checkbox"/>	Remote transmission alarm signal (via potential-free contact) in case of temperature alarm (change-over contact)
<input checked="" type="checkbox"/>	Designed and tested for climatic class SN (ambient temperature range +10°C to +32°C)
<input checked="" type="checkbox"/>	Interior made from stainless steel [evaporator plates and drip channel made from stainless steel]
<input checked="" type="checkbox"/>	Internal LED lighting
<input checked="" type="checkbox"/>	Additional remote transmission alarm signal (via potential-free contact) in case of system failure (change-over contact)
<input checked="" type="checkbox"/>	Smooth castors for optimum flexibility of movement
<input type="checkbox"/>	Water cooling [standard for water-cooled version: CSF61-101 W models]
<input type="checkbox"/>	Machine compartment as external unit [standard for split version: CSF61-101 S models]
<input checked="" type="checkbox"/>	Ambient temperature sensor
<input type="checkbox"/>	*B Connected - Universal software for the monitoring of refrigeration devices, including the acquisition, recording and visualization of temperature data
<input checked="" type="checkbox"/>	Ethernet interface for the visualization of all operating and control functions (hardware and software settings) via *B Connected monitoring software on a peripheral device (computer)
<input checked="" type="checkbox"/>	Battery operated emergency opening during the power failure

Contact Shock Freezers

B Medical Systems | CSF Range

2 models • Set temperature -30°C • Climate class SN • Compliant to 21CFR Part 864.9700, Class II | AABB Standards

In conformity with national and international guidelines, regulations for Medical Devices offering reliability, efficiency and safety at an optimal price.

MODELS CSF61-101 L

Fastest contact shock freezers worldwide

Our contact shock freezers provide the shortest freezing cycle on core temperature of -30°C : nearly 100 units of 250 ml within less than 60 minutes.

State-of-the-art contact freezing technology™

Our innovative and patent pending system of inclined trays delivers the most homogeneous freezing and offers the best in repeatability and productivity, with operator independent performance for 100% performance consistency.

Highest degree of ergonomics

Each tray has an effective and usable height, giving easier access for quicker loading and enabling easier manipulation of bags.



Optimized pull-down of the layers

The freezing down is quicker, the clinical payback on investment is maximised and the concentration of factors 8, 9 and albumin is increased.

Integrated multifunction electronics

Located at eye level, the electronics are easily accessible and can be connected with our fully integrated monitored solutions for a complete traceability of frozen bags.

Maintenance-friendly design

- Our contact shock freezers offer an easy cleaning and disinfection facility
- Thanks to the additional battery, plates can still be opened to access the bags, in case of power failure



Contact shock freezers are devices intended for the quick-freezing and subsequent temporary storage of blood plasma or biological samples to a core temperature below -30°C within less than one hour. CSF models are built in compliance with directives for the preparation of blood plasma and offers multiple controls, warning and data recording safety features developed as a result of 40 years of expertise in refrigeration technology.

CSF models are also available as split version (CSF61-101 S) and water-cooled version (CSF61-101 W). Technical data available upon request.



Technical Data

General features

CSF61 L


Freezing capacity (plasma bags)	<i>Single layer</i>	20 at 1000ml (content 850ml)
	<i>Double layer</i>	30 at 350ml (content 250ml)
Freezing time to core temperature of -30°C at +25°C ambient temperature (in single layer)		60 at 350ml (content 250ml)
		20 units (1000ml) ± 45 min.
Operating temperature (preset)	<i>Pre-cooling</i>	-40°C
	<i>Freezing</i>	-50°C
	<i>Eco-mode</i>	-43°C
Climate class (ambient temperature range)		SN (+10°C to +32°C)
Defrosting technique / time		Manual (hot gas) / 10 min.
Refrigerant type		R449a
External dimensions H x W x D (in)		75.2 x 42.5 x 32.7
Dimensions contact plates W x D (in)		2 units of 26.8 x 24.8
Net weight with standard equipment (lb)		1323

Energy Saver	Operating voltage range	380-400 V - 50/60 Hz (32A) or 480 V - 60 Hz (32A)		
	Power (W)	STABLE RUNNING unloaded	FREEZING PROCESS - fully loaded (with units at 1000 / 350ml in single layer)	
			8200	4800 / 3500 (5700 at defrosting)
	Energy consumption (kWh/24h)	0.7	3.0 / 1.8 (per freezing cycle)	
	Average heat ejection (Kcal/h)	-	3397	
	Compressor running time (%)	-	90	
	Noise level (dB(A)) (at 1.6m height & 1m distance)		73	73



OPTIMIZED CONTACT FREEZING

- Levels operate independently, allowing simultaneous freezing and defrosting as necessary
- Inclined plates at 5° ensure full contact with plasma and minimizes effect of air trapped in bags negatively influencing repeatability of freezing cycle times, providing operator independent performance for 100% performance consistency





Technical Data Specifications


CSF61 L

<i>Control/display sensor (4 sensors per device)</i> Precision (from -80°C to +180°C) ± 0.3°C	PT1000 2-WIRE 1/3DIN CL.B
<i>Dummy bag sensor (optional)</i> Precision (from -80°C to +180°C) ± 0.3°C	PT1000 2-WIRE 1/3DIN CL.B
<i>Accu data / function time of the control panel when system failure</i>	12 V - 7 AH / 48h
<i>Relative humidity at +32°C</i>	≤ 70%
<i>Compressor</i>	Bitzer
<i>Casing insulation (polyurethane)</i>	3.15 in PU
<i>Material</i>	
<i>Inner body</i>	Stainless steel (V2A - 1.4301)
<i>Outer body</i>	Steel (1.0330+ZE25/25) & Powder coated 60-90 μm
<i>Contact plates</i>	Stainless steel (V2A - 1.4301)
<i>European Medical Device Directive</i>	MDD 93/42/EEC, Class IIa
<i>EMC directive</i>	2014 / 30 / EU
<i>Low voltage directive</i>	2014 / 35 / EU
<i>FDA regulation Medical Device</i>	21CFR Part 864.9700, Class II

EXAMPLE OF "FREEZING CYCLE ON CORE TEMPERATURE OF -30°C"


CSF61 L

- Voltage: 400V / 50Hz
- Ambient temperature: +25°C
- Relative humidity: 45%
- Test bag type: 350ml, filled with 250ml (30 bags)
1000ml, filled with 850ml (20 bags)
- Filling medium: saline solution, 0.9% salt



CSF101 L

- Voltage: 400V / 50Hz
- Ambient temperature: +25°C
- Relative humidity: 45%
- Test bag type: 350ml, filled with 250ml (48 bags)
1000ml, filled with 850ml (32 bags)
- Filling medium: saline solution, 0.9% salt





DUMMY /
REFERENCE
BAG

(Shown with
sensor & cable)



TRANSPORT
CARRIER



TRANSPORT
TRAY



BARCODE
READER



EXTERNAL
AIR COOLING
SYSTEM

Equipment Standard & optional

CSF61 L

B Medical Systems Electronics 7" full touchscreen display	<input checked="" type="radio"/>	
Internal LED lighting	<input checked="" type="radio"/>	
GMP - Clean room classification	<input checked="" type="radio"/> A / ISO 5 <input type="radio"/> B / ISO 6	
°B Connected - Monitoring Software	<input type="radio"/>	
Ethernet interface / interface card	<input checked="" type="radio"/>	
Ambient temperature sensor	<input checked="" type="radio"/>	
Potential-free contact in case of system failure	<input checked="" type="radio"/>	
Eco-mode (interim storage temperature of -43°C)	<input checked="" type="radio"/>	
Dummy / reference bag port for accompanying core temperature readings during the freezing process	<input checked="" type="radio"/> 2	
Dummy / reference bag with sensor & cable	<input type="radio"/>	
Dummy / reference bag	<input type="radio"/> 500ml	
Transport carrier tray	<input type="radio"/> 1 <input type="radio"/> 4 (recommended)	
Barcode reader	<input type="radio"/>	
Condenser filter	<input checked="" type="radio"/>	
Integrated air-cooled condenser	<input checked="" type="radio"/>	
Water cooling	<input type="radio"/> <input checked="" type="radio"/> CSF61 W (water-cooled version)	
Air-cooled condenser & compressor as external unit (H 42.5 x W 53.1 x D 38.0 in)	<input type="radio"/> <input checked="" type="radio"/> CSF61 S (split version)	
Securing feet for models with condenser & compressor as external unit	<input type="radio"/>	
Smooth castors with stabilizers	<input checked="" type="radio"/>	
Wooden packaging for ocean transport / export	<input checked="" type="radio"/>	

INSTALLATION INSTRUCTIONS

The installation costs for units with an external air cooling system or water-cooled system are not included in the unit price (standard version: CSF61-101 L). These costs are subject to an individual client request and can vary depending on the conditions in the designated building. CSF models are installed as a free-standing, air-cooled unit ready for connection.

On request the CSF models can also be equipped with the following cooling systems:

- External air cooling system (CSF61-101 S models)
- Water cooling system (CSF61-101 W models)
 - Directly connected to the in-house water system
 - Connected to a separate closed air-conditioning system, also called "A/C water cooling system"
 - Connected to a separate closed water circuit, also called "refrigerated water tower"

Standard Optional - Not available

Full functionalities at a glance

B Medical Systems | Electronics

The new B Medical Systems Electronics

offers a wide range of adjustment and diagnostic facilities as well as additional protection / warning operations (via external alarm operations, histories and individual display signals).

Multi-function electronics with 7" full touchscreen display, located at eye level with pre-installed connection allowing exclusive °B Connected monitoring functionalities for a complete traceability of frozen bags.



Power switch OFF / ON

Turn main switch to ON position, the controller will be switched on automatically and the green LED lights up and the display activates.



Touchscreen 7" full display

The main screen has two different shapes depending on whether the batch management is active or not:

- At the top: login status (user level), current time and current voltage
- In the left-hand: the main screen shall always show two frames containing level specific data (operating mode, plate temperatures, dummy bag temperature, process time and progress bar). The entire frame and its content is an active touch area. When pressed, the selected level opens in detail view. The other level goes to minimized view
- In the right-hand: touch buttons to access batch list and main menu, to deactivate touch screen for cleaning, to mute alarm buzzer and to access active alarms and warnings. When the batch management is inactive, the batch list button is replaced by two moving direction selection button (open/close plates)



Modern and user-friendly graphic display

Acoustic and visual alarm signal

When an alarm occurs, the red alarm LED lights up and an acoustic signal sounds. The corresponding alarm message appears overall or level specific on the display. If several alarms occur simultaneously, the messages will be shown alternately.



Keycard reader

Login using a NFC card

Touch buttons

- CLEAN button: the cleaning button allows the user to disable the touch capability of the screen for a short period of time in order to clean the screen. This button is replaced by the logout button when any user is logged in
- MENU button: allows access to the menu screen
- MUTE button: deactivates the acoustic alarm for a predefined period of time
- ACTIVE ALARMS AND WARNING button: during a warning or alarm situation, the color of this button and the corresponding level will change. Information on an active alarm and warning is accessible by pressing the message button



SAVING LIVES THROUGH RELIABLE AND INNOVATIVE TECHNOLOGY

Optimum control and protection

B Medical Systems | °B Connected

°B CONNECTED | MONITORING SOFTWARE

Universal software for the monitoring of refrigeration devices, including the acquisition, recording and visualization of temperature data.



- Unique monitoring software for the full range of Blood Management Solutions and Medical Refrigeration products
- Web-based interfaces for computers and mobile devices
- Modern design for simple and intuitive use
- Graphic display of temperature curves
- Integrated event and activity history of appliances' components
- Data recording on centralized database for long-term archiving
- Easy setting of specific alarm, via email or SMS alerts
- Generation of reports compiling crucial data and activities
- Temperature and detailed device data export for third-party software
- Important cost advantage compared to a traditional circular chart recorder and its spare parts
- REST API to access raw data directly from the database in read-only mode

KEY BENEFITS:

- Unique interface for the monitoring of the full range of refrigeration devices
- Centralized database providing data access to entire customer network

1 CONFIGURATION




USER MANAGEMENT

- Password login
- Access right



DEVICE SETTINGS

- Device information
- Device & user configuration



ALARM PROFILE

- Parameters of the alarm profile (day, hour, minute)
- Alarm management

2 DISPLAYING



DETAILED DISPLAY

- Optimized design for mobile and touch compatible devices
- Support of multiple languages
 - Real-time data: Current temperatures, set point, upper and lower limit, door state, compressor runtime
- Log data: Events, alarms within date/time, event type, description of the event. Event and alarm confirmation function
- Batch management: Plasma bag barcodes scanned on device and imported in the database
 - Parameter settings of each device: Factory settings, real-time clock, language, firmware version, RS485 address, IP address

- °B Connected offer a highly flexible system that is adaptable to specific customer requirements
- > Complete & legally safe documentation of temperature data
- > Comprehensive applications and diagnostic possibilities

3 ANALYSING



REAL TIME GRAPHICS

- Graphic temperature progress (in real time)
- Adjustable time axis
- Offline navigation




INTERACTIVE ANALYSIS

- No limitation in terms of time range
- Freely scalable
- Detailed information for any point of time
- Temperature data, events and alarms
- Login history and audit trails



CONTROL CENTER

- Overview (Control Monitor) of all monitored appliances:
- Information regarding status, appliance name, location, connection, upper and lower limits, current temperature, events / alarms
- Dashboard to combine several temperature measurements in one single graph



PRINTED ANALYSIS

- Customizable print reports on daily, weekly and monthly basis
- Exports in various file types
- Standardized format