User Manual

Round Vacuum desiccators

Model: VDR-20(G/U/UG) and 25(G/U/UG) and VDR-30(G/U/UG)

Manual No.: D 2 011L002 Version : 3.0

Jeiotech

⚠ WARNING

Before using this product, read this entire Operator's Manual carefully. Users should follow all of the Operational Guidelines contained in this Manual and take all necessary safety precautions while using this product. Failure to follow these guidelines could result in potentially irreparable bodily harm and/or property damage.

Thank you for purchasing Jeio Tech's products.

Jeio Tech Co., Ltd. is committed to customer service both during and after the sale. If you have questions concerning the operation of your unit or the information in this manual, contact our Sales Department. If your unit fails to operate properly, or if you have questions concerning spare parts or Service Contracts, contact our Service Department.

Quality Management System

Jeio Tech Co, Ltd. is dedicated to providing world-best product quality and customer satisfaction. To ensure we maintain this commitment we have developed and implemented a total quality program, which conforms to the requirements according to DIN EN ISO 9001:2015 for the design, development, production, sales and servicing of biotechnology, environmental chemical engineering related products, and reliable measuring equipment for electric and electronics (ovens, incubators, constant temperature and Humidity chambers, constant temperature baths, refrigerating bath circulators, heat exchangers and shakers).

Visit our Web site at www.jeio.tech.com to view a copy of our certificate.

Disclaimer

Jeio Tech Co., Ltd. is committed to a continuing program of product development and improvement, and reserves the right to change information, such as specifications, appearance, and dimensions, described in this document without notice.

Copyright

No part of this manual may be reproduced or transmitted in any form or by any means, including photocopying, recording, or using information storage and retrieval systems, for any purpose other than the purchaser's own use, without the express written permission of Jeio Tech Co., Ltd.

©2025. All Rights Reserved. Jeio Tech Co., Ltd.

Any other product names and services identified in this manual are trademarks or registered trademarks of their respective owners. No such use, or the use of any trade name, is intended to convey endorsement or other affiliation with Jeio Tech Co., Ltd.

Table of Contents

1	Safety Advice	1
1.1	User Manual	1
1.2	DANGER/CAUTION/NOTICE Alerts.....	1
1.3	Danger Alerts	2
1.4	Caution Alerts	2
2	Product Description.....	3
2.1	Introduction	3
2.2	Characteristics.....	4
2.3	Layout.....	6
3	Unpacking and Installation	7
3.1	Checking Shipment Damage and Unpacking.....	7
3.2	Checking the Contents of the Package	7
3.3	Installation	8
3.3.1	Installation Environment.....	8
3.3.2	Installation location.....	8
4	Operation	9
4.1	Precautions.....	9
4.2	Operation of the 3-Way Valve.....	9
4.3	Vacuum formation	10
4.4	Storage	12
5	Maintenance.....	12
6	Troubleshooting.....	12
7	Accessories	13
8	Appendix	13
8.1	Technical Specifications	14
8.2	Disposal of the Unit	14
8.3	Warranty	15
8.3.1	Warranty period.....	15
8.3.2	Technical Service Contact Points.....	15
8.3.3	Certificate of Product Warranty	16

1 Safety Advice

1.1 User Manual

This manual contains important safety and operation information. You must carefully read, understand, and follow all the instructions in this manual prior to operating this instrument. Keep this manual in a safe place nearby for reference and make it easily available to all users.

1.2 DANGER/CAUTION/NOTICE Alerts

This manual highlights DANGER/CAUTION/NOTICE alerts to prevent injury or property damage and also to achieve optimum performance of your instrument.

These alerts are classified into three types in this manual depending on the importance and the risk levels as described below:

Symbols	Meaning
	Ignoring this alert could cause serious or even fatal injury.
	Ignoring this alert could cause serious injury or property damage.
	Ignoring this alert could cause operational problems.

Ensure to read the danger and caution alerts in 1.3 and 1.4 carefully before use.

1.3 Danger Alerts

 DANGER
Ignoring the following warnings could cause serious or even fatal injury.

- ◆ Never install or use this unit in explosive atmospheres.
- ◆ Never install this unit near to hazardous or flammable substances.
- ◆ Never disassemble, repair, or modify this unit on your own. Doing so will void your warranty and may result in injuries or product damages.

1.4 Caution Alerts

 CAUTION
Ignoring the following cautions could cause serious injury or property damage.

- ◆ Carefully read all the warning labels before using this unit. And do not remove or damage the warning labels.
- ◆ Do not move this unit during operation. When moving it, use the molded-in handles.
- ◆ Do not use this unit in wet or moist atmosphere or where water leakage is expected.
- ◆ Do not use this unit in contaminated atmosphere or where metallic dust exists.
- ◆ Do not expose this unit to any heat sources including direct sunlight.
- ◆ Do not allow moisture, organic solvent, dust, or corrosive gases to get into this unit.
- ◆ If you observe a strange smoke, odor or noise from this unit, stop using it. After the smoke or odor disappears, contact your dealer or Jeio Tech if any repair is required.

- ◆ Do not use chlorine bleach, ammonia-based cleaners, abrasives, ammonia, or metal scouring pads. Wipe with a soft damp cloth or a sponge soaked in water or diluted neutral detergent.

2 Product Description

2.1 Introduction

Congratulations of your purchase of our high-quality cylindrical vacuum desiccator which is specially designed to provide maximum benefit for your investment with respect to performance, safety, ease of use, and durability.

Currently eight cylindrical vacuum desiccator models of are available from Jeio Tech: **VDR-20(G/U/UG) and 25(G/U/UG) and VDR-30(G/U/UG)**

Note that the models with a **G** designation (i.e., VDR-20**G**/20**UG**/25**G**/25**UG**/30**G**/30**UG**) come standard with a vacuum gauge while the models with a **U** designation (i.e., VDR-20**U**/20**UG**/25**U**/25**UG**/30**U**/30**UG**) are UV blocking models. These UV blocking models are offered to minimize damages or discoloration of light-sensitive samples by using UV-blocking resin which completely blocks all ranges of UV and the blue spectrum of the visible light.

Being lightweight and yet shatter-resistant, Jeio Tech's Round vacuum desiccators are safer and economical alternatives to conventional glass desiccators for both vacuum and non-vacuum applications.

The vacuum environment created by these desiccators provides a very low humidity, contamination-free, corrosion-free, and dust-free condition which is ideal for:

- ◆ storing moisture-sensitive materials such as anhydrides and other hygroscopic compounds as well as materials which are highly reactive with atmosphere,
- ◆ removing residual solvents or air bubbles,
- ◆ cultivating anaerobic microbes,
- ◆ minimizing the reabsorption of humidity and the exposure to air during the drying or the cooling stages of biological/medical/food/materials experiments,
- ◆ conducting experiments at low pressure,
- ◆ drying at low pressure to minimize the heat effects.

2.2 Characteristics

(1) Outstanding Vacuum Sustainment

This unit maintains a 1 Torr (133 Pa) vacuum for over 72 hours at room temperature and stays gastight allowing reliable experiments over extended periods.

(2) Molded Polycarbonate Body and Base

The body and the base of this unit are made of shock-resistant, highly tensile, and transparent polycarbonate to provide extreme durability as well as convenient visual observation from outside.

The body is molded as a single piece requiring no adhesive and, therefore, allows virtually no leakage.

(3) Viton® 3-way Valve

This unit comes standard with a leak-proof Viton® 3-way valve offering consistent and uniform vacuum draw, vacuum release, or gas exchange without the inconvenience of connecting and disconnecting pump hoses to the unit.

If needed, an additional 3-way valve can be attached to the port at the base.

When releasing the vacuum, If using the Pore plug on the 3-Way Valve, it allows air to slowly enter the small diameter inlet, which protects dust or powder samples from flying out.

(4) Locker

This unit is equipped with a locker which allows tight sealing for both vacuum and non-vacuum applications. In addition, this locker provides additional convenience and safety when moving the unit.

(5) Wide Range of Models

Currently eight different models of are available: VDR-20(G/U/UG) and 25(G/U/UG) and VDR-30(G/U/UG)

These models come in two different capacities to suit the user's needs with the internal volume of either 10 liter (for VDR-25 models) or 20liter (for VDR-30 models).

The models with a **G** designation (i.e., VDR-20**G**/20**UG**/25**G**/25**UG**/30**G**/30**UG**) come standard with a vacuum gauge while the models with a **U** designation (i.e., VDR-20**U**/20**UG**/25**U**/25**UG**/30**U**/30**UG**) are UV blocking models.

(6) Easy Vacuum Level Checking

With a vacuum gauge installed on top of the body, the vacuum level checking can be done conveniently for the models equipped with a vacuum gauge (VDR-25G/25UG/30G/30UG).


(7) UV Blocking

These UV blocking models (VDR-20U/20UG/25U/25UG/30U/30UG) are offered to minimize damages or discoloration of light-sensitive samples by using UV-blocking resin which completely blocks all ranges of UV and the blue spectrum of the visible light.

(8) High Quality Vacuum Seal

The silicone vacuum seal ensures reliable vacuum formation and sustainment without using any grease. For easy cleaning and maintenance, the vacuum seal is specially designed to fit snugly into the bottom groove of the body.

2.3 Layout

[VDR-30G Model]

1. Desiccator jar
2. Flange Locker
3. Vacuum Seal
4. Desiccator cover
5. Vacuum gauge
6. 3-Way Valve
7. Perforated sample tray
8. Desiccant tray


3 Unpacking and Installation

3.1 Checking Shipment Damage and Unpacking

Upon receiving the instrument, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred during shipment is detected before unpacking. If you find such damage, the carrier must be notified immediately.

3.2 Checking the Contents of the Package

After unpacking, check to ensure that all the parts and accessories described below are included in the package. If not, contact your dealer immediately.

Item	Figure	Quantity	Description
Main Body	-	1	-
3-way valve		1	-
Vacuum Gauge (Option)		1	Only G series
Pore plug		1	-
Vacuum seal		1	-
Locker		1	-
Perforated Shelf	-	1	-
Desiccant tray	-	1	-

User manual		1	-
-------------	---	---	---

3.3 Installation

3.3.1 Installation Environment

NOTICE
Avoid exposure to any heat sources including direct sunlight.

This unit is designed for indoor use in laboratory environments. The conditions for proper installation are as follows:

Permissible Ambient temperature	5 to 40°C (41 to 104° F)
Permissible Relative humidity	below 80%
Permissible Altitude	0 to 2,000 m (6,562 ft)

3.3.2 Installation location


NOTICE
Install the unit on a level surface.

4 Operation

4.1 Precautions

 CAUTION
• Do not form a vacuum when airtight containers are inside the unit. • Note that this unit cannot be autoclaved. • When attaching other equipments or additional 3-Way Valves after removing the port plugs, use proper O-rings and nuts for tight sealing. • When putting high-temperature samples inside the unit, cool them thoroughly in advance or make sure that they are not in direct contact with the interior of the unit to avoid interior damage.

4.2 Operation of the 3-Way Valve


ITEM	Front Nozzle	Rear Nozzle
	Open	Closed
	Closed	Open
	Closed	Closed

4.3 Vacuum formation

The procedures for vacuum formation are as follows:

- (1) Connect the vacuum pump hose to one of the two nozzles of the 3-Way Valve.
- (2) After activating the vacuum pump, turn the valve handle towards the nozzle to open it and read the vacuum gauge to check whether the vacuum is forming properly.
- (3) When the desired vacuum level is achieved, turn the valve handle until it is perpendicular to the nozzle to close it and check the vacuum gauge whether there is any leakage.
- (4) When releasing the vacuum, slowly turn the valve handle towards the opposite nozzle to return to the normal pressure.
→ At this time, using the Pore plug on the 3-Way Valve allows air to slowly enter the small diameter ($\varnothing$ 1 mm) inlet, which protects dust or powder samples from flying out.

Here's how to use the Pore plug on the 3-Way Valve.

The cap of Pore plug

- ① Insert the loop, the end part of the Pore plug into the nozzle(Air inlet) of 3-Way Valve.
- ② Cap the Pore plug into the nozzle(Air inlet) of 3-Way Valve.
- ③ Turn the valve towards the capped air inlet to release the vacuum.

However, if you use the Pore plug, it takes 3 to 5 times longer than using 3-Way Valve only to recover to atmospheric pressure. The following table compares the recovery time

to atmospheric pressure when using the Pore plug on 3-Way Valve for each model. Please, use it with reference.

Model		VDR-20G	VDR-25G	VDR-30G
Internal volume		6 L	10 L	20 L
Time to recover to atmosphere pressure	3-Way Valve Air inlet (Ø, 4.5mm)	About 10~ 30sec	About 30 sec	About 40 sec ~ 1 min
	Pore plug (Ø, 1mm)	About 1~2 min	About 2~3 min	About 4~5 min

CAUTION

- Do not impact or drop the unit when a vacuum is formed inside.
- When forming or releasing a vacuum, turn the valve handle slowly to avoid rapid pressure change. If not, the stored samples can be damaged or dispersed.
- When forming a vacuum, install a cold trap in front of the vacuum pump. Moisture and chemicals from the samples may result in damages to the pump.
- Do not let the body of the unit into direct contact with acetone, benzene, toluene, chloroform, cresol, sodium hydroxide, highly concentrated nitric or sulfuric acid, acetic acids, or strong chlorine-based solvents.
- Ensure that the door seal is not damaged by strong acid or base. Damaged seal can cause malfunction of the unit.
- When a vacuum is formed inside the unit, do not apply excessive force to the 3-Way Valves, the vacuum gauge, or any of their vicinities. Any damages caused by external forces can cause malfunction of the unit.
- When the unit is not forming a vacuum, keep the valve handle perpendicular to the nozzles.
- When relocating the unit, do not hold on to protruded parts such as the vacuum gauge or 3-way valves.

4.4 Storage

If this unit is not used for an extended period of time, clean the instrument with soft cloth and store it in dry place after packing it properly.

5 Maintenance

 CAUTION	
● If the unit is contaminated, wear chemicals-resistant gloves before cleaning. ● Regularly check the 3-Way Valves, the vacuum gauge, or any of their vicinities for any damage. ● Regularly check the door seal as well as the surface where the door and the body are in contact and keep them always clean and undamaged. ● Make sure to clean the door seal regularly using diluted neutral detergent. For thorough cleaning, detach the seal from the door. ● Do not use chlorine bleach, ammonia-based cleaners, abrasives, organic solvents, or metal scouring pads when cleaning. Use a soft cloth all the time.	

6 Troubleshooting

Problem	Cause	Corrective Action
Vacuum Malfunction	Contaminated surface between the door and the body	Thoroughly clean the surface with a soft cloth using diluted neutral detergent.
	Damaged vacuum seal	Replace the seal.
	Unsuitable room temperature	Check the temperature whether it is within the operating range.

	Misplaced or damaged 3-Way Valve	Check whether the valve is properly set. If damaged, replace it with a new one.
	Damaged vacuum gauge	Replace the gauge.

7 Accessories

Component	Cat. No.	Description	
Vacuum Gauge	EDA4300	0 -0.1(Mpa)	
Desiccant Tray	AAAD150 1	VDR-20(G,U,UG)	-
	AAAD1502	VDR-25 (G,U,UG)	
	AAAD1503	VDR-30 (G,U,UG)	
Perforated shelf	AAAD1511	VDR-20(G,U,UG)	-
	AAAD1512	VDR-25 (G,U,UG)	
	AAAD1513	VDR-30 (G,U,UG)	
Vacuum seal	VDR0001	VDR-20(G,U,UG)	Silicone Rubber
	VDR0002	VDR-25 (G,U,UG)	Silicone Rubber
	VDR0003	VDR-30 (G,U,UG)	Silicone Rubber
Locker	VDR0015	VDR-20(G,U,UG)	-
	VDR0016	VDR-25 (G,U,UG)	-
	VDR0017	VDR-30 (G,U,UG)	-
3-Way Valve	AAAD1521	PP, Ø9.5 x 18	
Silica Gel	BEA00078603	Silica Gel (1kg)	

8 Appendix

8.1 Technical Specifications

MODEL		VDR-20(G)	VDR-25(G)	VDR-30(G)
		VDR-20U(UG)	VDR-25 U(UG)	VDR-30 U(UG)
Internal volume (L)		6	10	20
Vacuum range (MPa)		0 ~ -0.1 (Analog)		
Nozzle (mm)		Ø 9.5		
Max. permissible vacuum (Bodies)		1.33 x 10 ⁻⁴ MPa		
Material	Cover/Jar	PC (Polycarbonate)		
	Tray, Locker	PP (Polypropylene)		
	Vacuum seal	Silicone rubber		
Dimensions (mm, WxH)	Overall(mm, WxH)	286 x 279 (286x 354)	353 x 325 (353 x 400)	427 x 399 (427x 475)
	Net weight(kg)	1.5(1.7)	2.5(2.6)	4.2(4.3)
Shelves	Perforated Shelf(mm)	192	250	310
	Max load(Kg)	15	20	25

Model	VDR-20U / 20UG	VDR-25U / 25UG	VDR-30U / 30UG
Color	Amber		
UV Protection (Range)	UV A,B,C 100% (200~450nm)		

8.2 Disposal of the Unit

Disposing of the unit must be done in an environmentally responsible way if it has been

potentially exposed to bio-agents or radioactive samples. Failure to follow stringent requirements for appropriate disposal may lead to actions against you and your organization.

- (1) First, check with your laboratory or organization to ensure that you are following all the policies and procedures for disposal of laboratory equipments and wastes.
- (2) If not possible, contact your local governing body for regulations regarding disposal of laboratory equipments. Jeio Tech highly recommends you to find a local service provider that can properly dispose of your unit.

8.3 Warranty

8.3.1 Warranty period

- (1) If the product fails during normal and proper use, the warranty period for manufacturing liability shall be two years from the date of delivery.
- (2) If you request a repair, please check the items below and let us know so that you can get a quicker and more accurate repair.

- Date of purchase:
- Customer name / address / phone number / E-mail:
- Fault status
- Model
- Serial number

8.3.2 Technical Service Contact Points

Technical Services

Address : 153 (Youngsan-dong), Techno 2-ro Yuseong-gu, Daejeon-city, 34025, Republic of Korea

Tel : 1522-4298

Website : www.jeiotech.com

E-mail : help@jeiotech.com

8.3.3 Certificate of Product Warranty

This is a guide to the product quality warranty and service terms provided by JEIOTECH Co., Ltd

All products are guaranteed to operate normally when installed under proper conditions, as specified in the user manual for each product, and used according to the intended purposes and permissible operating conditions.

The warranty period for product quality starts from the delivery date and is indicated in the catalog provided for each product or in the purchase order agreed upon with the customer. If a product malfunctions within the warranty period, JEIOTECH Co., Ltd. provides free repair services. However, among the transportation costs incurred during the repair process, any special costs resulting from the customer's site conditions (e.g., door or wall modifications, use of special handling equipment) are not covered under the warranty. Furthermore, JEIOTECH's liability under this warranty is limited to replacing the product or refunding the original purchase price at JEIOTECH's sole discretion if reasonable attempts to repair or replace the product fail.

Parts replaced during free repairs within the warranty period become the property of JEIOTECH Co., Ltd. Additionally, the warranty period is not extended due to free repairs. JEIOTECH is not responsible for any losses incurred by the customer due to delays in repairs caused by issues such as part shortages or shipping delays.

The following cases will be handled as paid repairs, even within the warranty period:

1. Non-defective cases

- * Requests to resolve issues arising from improper operating environments or abnormal installations.
- * Reinstallation requests or connections with external equipment other than the initial installation.
- * Requests for usage instructions or simple adjustments that do not require disassembly.
- * Inspections of non-defective products.
- * Cleaning or debris removal service requests.

2. Failures caused by user errors:

- * Malfunctions due to usage not in accordance with the user manual.
- * Issues resulting from improper operating environments or abnormal installations.
- * Malfunctions caused by incompatible chemical usage.
- * Malfunctions caused by irregular or improper electrical or utility supplies.
- * Damage due to product movement, storage, external impact, or dropping by the customer.
- * Malfunctions resulting from unauthorized modifications or repairs without prior written approval from JEIOTECH Co., Ltd.
- * Malfunctions caused by using consumables or optional components not specified by JEIOTECH Co., Ltd.
- * Repairs performed by personnel other than JEIOTECH or its authorized service agents.

3. Other cases

- * When the product model or serial number has been damaged, removed, manipulated, or altered.
- * Malfunctions caused by fire, flooding, or other natural disasters.
- * Failure due to consumable parts reaching the end of their lifespan (e.g., lamps, filters, fuses,

gaskets, packings, hoses, heat transfer fluids).

* Natural wear and tear (e.g., discoloration, screen printing degradation).

* Damage or malfunction of product accessories (e.g., power cables, data cables, glass components) or options.

* Malfunctions caused by errors in computers, computer components, or connected network systems.

- Before JEIOTECH performs free or paid repairs or accepts the return of a product, the customer must ensure the product has been decontaminated of chemical and biological hazards.
- Replacement parts used during repairs may differ from the original but will function correctly and be equivalent in performance.
- JEIOTECH does not bear any responsibility for damage to stored samples, data loss, or any other direct or indirect damages caused by product malfunctions, regardless of the warranty period.
- JEIOTECH does not guarantee issues related to computers, computer components, or network systems connected to the product, regardless of the warranty period.
- This warranty is the sole and exclusive warranty provided. Any implied warranties imposed by applicable laws are limited to the warranty period stated herein. This warranty cannot be modified by any verbal advice, statements, or documents on behalf of JEIOTECH unless explicitly labeled as a "Warranty Amendment" and signed by an authorized representative of JEIOTECH.
- This warranty applies only to products purchased for direct use by the customer. It does not cover products purchased for resale purposes without prior written approval from JEIOTECH.
- This warranty does not apply to consumables and is non-transferable.

- ✧ The contents of this manual can be changed or upgraded without prior notice.
- ✧ The copyright of this manual is reserved by Jeio Tech.

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