

Freezers

Undercounter and Counter-top

Marvel Scientific freezers are designed and manufactured to consistently meet the needs of health, science and industry. Marvel Scientific freezers have the power to cool contents quickly and dependably. Select freezers ranging in size from the 1.5 cu.ft. counter-top freezer to the large 6.1 cu.ft. undercounter model.

1 15AF General Purpose

- -20° C counter-top freezer
- Flat door liner
- Compact footprint
- Counter-top point-of-use storage or built-in undercounter

2 4CAF General Purpose

4FAF Flammable Material 

4EAF Hazardous Location 

- Flat door liner
- Removable shelf

3 6CAF General Purpose

- 3 door shelves
- 1 adjustable shelf
- 1 fixed shelf



Model Number	Volume cu.ft./liters	Interior Dimensions H x W x D in/cm	Exterior Dimensions H x W x D in/cm	Temperature Range °F/°C	Defrost	Shipping Weight lbs/kg
--------------	-------------------------	--	---	-------------------------------	---------	---------------------------

General Purpose Refrigerators

1 15AF	1.5 cu.ft. 42.5 liters	12 x 11¼ x 17½ 30 x 29 x 45	24¼ x 14⅞ x 22½ 62 x 38 x 56	-4° -20°	manual	70 32
2 4CAF	4.5 cu.ft. 127 liters	24¾ x 19 x 18¼ 63 x 48 x 46	34 x 23⅞ x 24 86 x 61 x 61	-10° to 0° -23° to -18°	manual	120 55
3 6CAF	6.1 cu.ft. 173 liters	27½ x 21 x 21 70 x 53 x 53	34 x 23⅞ x 24 86 x 61 x 61	10° to 20° -12° to -7°	manual	120 55

Flammable Material

2 4FAF*	4.5 cu.ft. 127 liters	24¾ x 19 x 18¼ 63 x 48 x 46	34 x 23⅞ x 24 86 x 61 x 61	0° to 5° -18° to -15°	manual	120 55
---------	--------------------------	--------------------------------	-------------------------------	--------------------------	--------	-----------

Hazardous Location

2 4EAF*	4.5 cu.ft. 127 liters	24¾ x 19 x 18¼ 63 x 48 x 46	34 x 23⅞ x 24 86 x 61 x 61	0° to 5° -18° to -15°	manual	120 55
---------	--------------------------	--------------------------------	-------------------------------	--------------------------	--------	-----------

*Note: Flammable material storage and hazardous location models feature static condensers and need an additional 1" clearance from top and sides and 1½" from back for proper condenser ventilation when recessed.