

# Laboratory Refrigerated Incubators


Model BOD05



Model BOD05



- Sizes range from 5.0 to 20.7 cubic feet
- Variable temperature range
- Microprocessor standard for better temperature control and stability
- Uniformity better than +/- .75°C



| Catalog Number | Type         | Cubic Feet | Temperature °C | Defrost | Doors   | Shelves      | Interior Dimensions W" D" H" | Exterior Dimensions W" D" H" | Shipping Weight (lbs.) | Commercial UL | Energy Star | Temperature Controls | Hex Lock |
|----------------|--------------|------------|----------------|---------|---------|--------------|------------------------------|------------------------------|------------------------|---------------|-------------|----------------------|----------|
| BOD05          | Undercounter | 5.0        | 0 to 50        | Manual  | 1 swing | 2 adjustable | 18.6 18.5 27.6               | 21.7 24.7 33.5               | 99.0                   | -             | -           | Microprocessor       | -        |
| BOD21          | Upright      | 20.7       | -9.9 to 50     | Manual  | 1 swing | 4 fixed      | 27.6 22.8 56.8               | 33.0 30.5 66.0               | 211.0                  | -             | -           | Microprocessor       | *        |