

# Green Innovation

## LABREPCO Low Temp Freezers (-86°C)

	<b>LABREPCO Harris Plus</b>	<b>Leading Competitor</b>
<b>Annual Energy Consumption</b>	<b>7,775 kWh</b>	<b>9,308 kWh</b>
<b>Annual Operating Cost Comparison</b>	<b>\$709</b>	<b>\$849</b>
<b>Life Cycle Energy Costs</b>	<b>\$7,090</b>	<b>\$8,490</b>

**Requires 16% Less Energy to Operate**

### Other Savings:

15% less heat rejection into room  
save on HVAC costs to cool room

CFC-Free Refrigerants and Insulation

40% quieter than previous designs



View our full line of Low Temp Freezers



View Our Other Eco-Friendly Products

**LABREPCO Harris PLUS  
Low Temp (-86°C) Freezer**

