

MetroMax® Lab Worktables — 52.01

- Advanced polymer construction provides excellent corrosion-resistance performance for the most corrosive lab environments.
- Table features smooth surfaces that are easy-to-clean. Mobile units can be easily moved to clean behind table area.
- Available top materials include Type 304 stainless and TRESPA phenolic resin (black or gray). Refer to chart below for more information.
- Choose from a wide range of available accessories and caster options (page 192) to customize table for a specific application.

Countertop Application Matrix							
Use Area	Epoxy Resin	Phenolic Resin	Stainless Steel	Characteristics	Epoxy Resin	Phenolic Resin	Stainless Steel
	2			Heat Resistance	Excellent	Fair*	Excellent
Low chemical use				Abrasion Resistance	Good	Excellent	Good
Moderate chemical use	X	X		Impact Resistance	Good	Excellent	Good
High chemical use		X	X	Strength	Good	Excellent	Good
				Stain Resistance	Good	Good	Good

*Discolors with high heat

Note: Phenolic resin and stainless are offered as top options for this table system. **Epoxy resin information is shown for comparison only.**



LTM60XUPG
(Shown with accessories, casters and StarSys Cart)

Definitions:

Phenolic resin: Compression molded of organic fiber-reinforced phenolic core with an integrally cured thermoset resin surface. Stronger than Epoxy resin over a span, and not as brittle.

Stainless Steel: Worksurfaces are fabricated from Type 304 stainless steel which provides excellent protection in corrosive environments.

Epoxy resin: Epoxy resins, silica, organic fillers, and inert hardeners are cast in forms and cured to a uniform mixture throughout.

Select Table Style from Chart Below and Following Page:

Worktables with Stainless Top and Solid MetroMax Shelf

Width (in.) (mm)	Length (in.) (mm)	Work Surface Height (in.) (mm)	Overall Height (in.) (mm)	Approx. Pkd. Wt. (lbs.) (kg)	Cat. No. Stationary	Cat. No. Mobile*		List Price Each				
30	762	29 ³ / ₄	749	35 ¹ / ₂	902	39 ¹ / ₂	1003	65	29	LT30XS	LTM30XS	962.00
30	762	35 ³ / ₄	908	35 ¹ / ₂	902	39 ¹ / ₂	1003	70	32	LT36XS	LTM36XS	1,005.50
30	762	47 ³ / ₄	1213	35 ¹ / ₂	902	39 ¹ / ₂	1003	90	41	LT48XS	LTM48XS	1,059.50
30	762	59 ³ / ₄	1517	35 ¹ / ₂	902	39 ¹ / ₂	1003	100	45	LT60XS	LTM60XS	1,102.50

*Casters ordered separately. Refer to page 195 for caster options.

Note: Lab Worktables are load rated at 50 lbs. (23kg) per square foot up to a maximum of 600 lbs. (273kg), assuming load evenly distributed.



LTM30XS
(Featuring solid MetroMax shelf shown with casters)

Worktables with Stainless Top and 3-Sided Frame

Width (in.) (mm)	Length (in.) (mm)	Work Surface Height (in.) (mm)	Overall Height (in.) (mm)	Approx. Pkd. Wt. (lbs.) (kg)	Cat. No. Stationary	Cat. No. Mobile*		List Price Each				
30	762	29 ³ / ₄	749	35 ¹ / ₂	902	39 ¹ / ₂	1003	70	32	LT30XUS	LTM30XUS	865.00
30	762	35 ³ / ₄	908	35 ¹ / ₂	902	39 ¹ / ₂	1003	75	34	LT36XUS	LTM36XUS	886.50
30	762	47 ³ / ₄	1213	35 ¹ / ₂	902	39 ¹ / ₂	1003	95	43	LT48XUS	LTM48XUS	919.00
30	762	59 ³ / ₄	1517	35 ¹ / ₂	902	39 ¹ / ₂	1003	110	50	LT60XUS	LTM60XUS	962.00

*Casters ordered separately. Refer to page 195 for caster options.

Note: Lab Worktables are load rated at 50 lbs. (23kg) per square foot up to a maximum of 600 lbs. (273kg), assuming load evenly distributed.



LTM30XUS
(Shown with accessory wire shelf and casters)

MetroMax Lab Worktables — 52.01

Select Table Style from Chart Below:

Worktables with **Black Phenolic Top and Solid MetroMax Shelf**

Width (in.)	Length (mm)	Work Surface Height (in.)	Overall Height (in.)	Approx. Pkd. Wt. (lbs.)	Cat. No. Stationary	Cat. No. Mobile*	List Price Each
30	762	29 ³ / ₄ 749	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	85 38	LT30XPB	LTM30XPB 940.50
30	762	35 ³ / ₄ 908	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	91 41	LT36XPB	LTM36XPB 983.50
30	762	47 ³ / ₄ 1213	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	118 53	LT48XPB	LTM48XPB 1,037.50
30	762	59 ³ / ₄ 1517	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	138 62	LT60XPB	LTM60XPB 1,081.00

*Casters ordered separately. Refer to chart below for caster options.
Note: Lab Worktables are load rated at 50 lbs. (23kg) per square foot up to a maximum of 600 lbs. (273kg), assuming load evenly distributed.



LTM60XPB
(Shown with back splash, accessories and casters)

Worktables with **Black Phenolic Top and 3-Sided Frame**

Width (in.)	Length (mm)	Work Surface Height (in.)	Overall Height (in.)	Approx. Pkd. Wt. (lbs.)	Cat. No. Stationary	Cat. No. Mobile*	List Price Each
30	762	29 ³ / ₄ 749	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	90 41	LT30XUPB	LTM30XUPB 843.00
30	762	35 ³ / ₄ 908	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	97 44	LT36XUPB	LTM36XUPB 865.00
30	762	47 ³ / ₄ 1213	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	125 56	LT48XUPB	LTM48XUPB 897.00
30	762	59 ³ / ₄ 1517	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	145 66	LT60XUPB	LTM60XUPB 930.00

*Casters ordered separately. Refer to chart below for caster options.
Note: Lab Worktables are load rated at 50 lbs. (23kg) per square foot up to a maximum of 600 lbs. (273kg), assuming load evenly distributed.



LTM60XUPB
(Shown with accessories and casters)

Worktables with **Gray Phenolic Top and Solid MetroMax Shelf**

Width (in.)	Length (mm)	Work Surface Height (in.)	Overall Height (in.)	Approx. Pkd. Wt. (lbs.)	Cat. No. Stationary	Cat. No. Mobile*	List Price Each
30	762	29 ³ / ₄ 749	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	85 38	LT30XPG	LTM30XPG 940.50
30	762	35 ³ / ₄ 908	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	91 41	LT36XPG	LTM36XPG 983.50
30	762	47 ³ / ₄ 1213	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	118 53	LT48XPG	LTM48XPG 1,037.50
30	762	59 ³ / ₄ 1517	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	138 62	LT60XPG	LTM60XPG 1,081.00

*Casters ordered separately. Refer to chart below for caster options.
Note: Lab Worktables are load rated at 50 lbs. (23kg) per square foot up to a maximum of 600 lbs. (273kg), assuming load evenly distributed.



LTM60XPG
(Shown with back splash, accessories and casters)

Worktables with **Gray Phenolic Top and 3-Sided Frame**

Width (in.)	Length (mm)	Work Surface Height (in.)	Overall Height (in.)	Approx. Pkd. Wt. (lbs.)	Cat. No. Stationary	Cat. No. Mobile*	List Price Each
30	762	29 ³ / ₄ 749	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	90 41	LT30XUPG	LTM30XUPG 843.00
30	762	35 ³ / ₄ 908	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	97 44	LT36XUPG	LTM36XUPG 865.00
30	762	47 ³ / ₄ 1213	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	125 56	LT48XUPG	LTM48XUPG 897.00
30	762	59 ³ / ₄ 1517	36 ¹ / ₂ 927	40 ¹ / ₂ 1016	145 66	LT60XUPG	LTM60XUPG 930.00

*Casters ordered separately. Refer to chart below for caster options.
Note: Lab Worktables are load rated at 50 lbs. (23kg) per square foot up to a maximum of 600 lbs. (273kg), assuming load evenly distributed.



LTM60XUPG
(Shown with accessories and casters)

For Mobile Worktables, Select Appropriate Casters from Chart Below. (Four brake casters are recommended for maximum table stability)

Caster Selection Guide

*Caster Series	Wheel Material	Rollability	Floor Protection	Noise	Load Rating (lbs.)	Load Rating (kg)	Corrosion Resistance	Cart Washable	Antimicrobial Additive in Wheel
5MDXA/5MDBXA (brake)	High Modulus Rubber	Good	Good	Low	250	114	Low	No	No
5MPX/5MPBX (brake)	Polyurethane	Good	Good	Moderate	300	137	Low	No	No
5PCX/5PCBX (brake)	Polyurethane	Good	Good	Moderate	300	137	High	No	No
5PCXM/5PCBXM (brake)	Polyurethane	Good	Good	Moderate	300	137	High	No	Yes
5MPXGSA/5MPXGSA (brake)	Polyurethane	Good	Good	Moderate	300	137	High	Yes	No

NOTE: The total weight of the equipment and its load should not exceed three times the load rating per caster. Additional stem casters, in various sizes and wheel materials, are available. Contact your Metro rep for more information.



5MDXA



5MPX



5PCBX



5PCBXM



5MPXGSA

Note: Use heavy-duty mobile post adapter (Cat. No. 9974HDP-4) to convert HD Mobile post to accept any standard Metro 5" (127mm) series casters (shown above).

MetroMax Lab Worktable Accessories — 52.01

Backsplash

4" (102mm) high. Attaches to phenolic table top (page 195), using pre-drilled holes. Match length and color with selected countertop.

Length (in.) (mm)		Approx. Pkd. Wt. (lbs.) (kg)		Cat. No. Black	Cat. No. Gray	List Price Each
27 ³ / ₄	705	5.5	2.5	LT30BSB	LT30BSG	65.00
33 ³ / ₄	875	6.5	2.9	LT36BSB	LT36BSG	78.00
45 ³ / ₄	1162	7.5	3.4	LT48BSB	LT48BSG	95.50
57 ³ / ₄	1466	9.5	4.3	LT60BSB	LT60BSG	113.00



L30BSG

Keyboard Tray

Fits 48" and 60" (1219 and 1524mm) table lengths. Cat. No. **LTKB** List Price 235.00 each

Note: Keyboard tray and drawer will both fit a 60" (1524mm) table.



Drawers

Fit 30", 36", 48" and 60" (760, 914, 1219 and 1524mm) table lengths.

3" (76mm)	Cat. No. LTD3	List Price 195.00 each
6" (152mm)	Cat. No. LTD6	List Price 255.50 each
3" (76mm) locking	Cat. No. LTD3L	List Price 281.50 each
6" (152mm) locking	Cat. No. LTD6L	List Price 342.00 each

Note: Keyboard tray and drawer will both fit a 60" (1524mm) table.



Wire Shelf

Stainless or Epoxy-coated shelves available.

Length (in.) (mm)		Approx. Pkd. Wt. (lbs.) (kg)		Cat. No. Stainless	List Price Each	Cat. No. Epoxy Coated	List Price Each
22 ³ / ₈	568	6	2.7	LT30WS	164.00	LT30WE	72.00
28 ³ / ₈	720	7.5	3.4	LT36WS	186.00	LT36WE	85.00
40 ³ / ₈	1025	10	4.5	LT48WS	207.50	LT48WE	98.00
52 ³ / ₈	1330	12	5.4	LT60WS	229.50	LT60WE	108.50

Note: Shelves are load rated for a maximum of 150 lbs. (68kg).



Surface Mount Swing Arm for Flat Monitor

Cat. No. **LTFMA**
List Price 487.00 each



LTFMA



MetroMax® Lab Carts — 52.01

Metro's Lab Cart line includes a General Lab Cart, Chemical Cart, and Mobile Desk. All carts are designed to offer superior corrosion resistance. Smooth rolling, 5" (127mm) polyurethane casters come standard with each cart. All shelves adjust in 1/2" (12.7mm) increments, offering precise shelf adjustment and maximum storage utilization.



LABX2
General Lab Cart
(casters included)



MDSX2
Mobile Desk
(casters included)



CHEMX2
Chemical Cart
(casters included)

Width (in.)	(mm)	Length (in.)	(mm)	Height (in.)	(mm)	Cart Load Rating		Description	Approx. Pkd. Wt.		Cat. No.	List Price Each
						(lbs.)	(kg)		(lbs.)	(kg)		
18	457	36	914	40	1016	600	273	General Lab Cart	43	19.5	LABX2	628.00
18	457	36	914	40	1016	600	273	Chemical Cart	49	22.2	CHEMX2	712.50
18	457	36	914	40	1016	600	273	Mobile Desk	55	24.9	MDSX2 ¹	1,168.50

¹ Mobile desk features a phenolic work surface offering a high level of corrosion resistance. Mobile desk includes: four polyurethane brake casters, each measuring 5" (127mm), for maximum stability when foot-operated brakes are engaged.
Note: All carts include 5MPX series casters (2 brake, 2 swivel), unless otherwise noted.

Deep Ledge Lab Solutions — 52.01

(with *Microban® Antimicrobial Product Protection)

Metro Deep Ledge Carts with Microban® antimicrobial product protection are the most advanced general purpose carts in the marketplace today. The Deep Ledge series has shelves with built-in *Microban® antimicrobial product protection that helps keep shelves cleaner between cleanings by inhibiting the growth of bacteria, mold and mildew that cause odors and stains on the shelf surface.

All models listed below come without casters. Refer to page 26 to select casters.

Width (in.)	(mm)	Length (in.)	(mm)	Height (in.)	(mm)	Description	Approx. Pkd. Wt.		Cat. No.	List Price Each
							(lbs.)	(kg)		
21 1/2	546	32 3/4	831	41	1041	Lab Dolly with Stainless Posts and Handle	25	11	LDP	410.00
21 1/2	546	32 3/4	831	41	1041	2-Shelf Cart with Stainless Posts and Handle	29	13	LC2P2030	323.50
21 1/2	546	32 3/4	831	41	1041	3-Shelf Cart with Stainless Posts and Handle	38	17	LC3P2030	432.00

Note: Cat. No. LC3P2030 is load rated for a maximum of 200 lbs. (91kg).



LC3P2030
(Shown with stainless casters)



LDP
(Shown with stainless casters)

*MICROBAN® and the MICROBAN® symbol are registered trademarks of the Microban Products Company, Huntersville, NC.