

STAINLESS STEEL MOBILE TRANSPORT CAGES



102305 Mobile Oxygen Therapy Unit

This Stainless Steel Mobile Transport Cage includes an oxygen therapy door. The door is made of 1/4" thick clear acrylic and stainless steel hardware. The door is complete with continuous flow nebulizer, ice compartment, thermometer and hydrometer.

Overall Size: 36"W 36"H 27.68"D



107373 Mobile Transport Cage

This Stainless Steel Mobile Transport Cage has an opaque door that is made of dark, tinted, 1/4" thick, acrylic, with ventilation holes. The door frame and hardware are made of stainless steel.

Overall Size: 36"W 36"H 25.5"D



102295 Mobile Transport Cage

This Stainless Steel Mobile Transport Cage has a stainless steel rod wire door which is electropolished.

Overall Size: 36"W 36"H 25.5"D

Each of these three mobile transport cages has a stainless steel body, featuring heli-arc welded construction with a .750" radius for all inside bends and coved corners.

Each of these cages has a flat stainless steel top, inside panels, and push handle. They are also offered with three different types of doors.

The cage is equipped with chrome plated swivel casters with 5" polyolefin wheels. Two of them have brakes. All wheels have swivel and axle bearings with seals and zerks grease fittings.

SSCI reserves the right to change any dimensions.

wwwsuburg.com

OPTIONS

Stainless Steel casters
Plastisol coated expanded metal floor
Stainless Steel floor
Stainless Steel collection pan
Stainless Steel shelf



STAINLESS STEEL MOBILE TRANSPORT CAGES with REMOVABLE COLLECTION PANS



Mobile Transport Cage with two rod grill hinged lids with rod grills welded in the front, sides, and back with optional bumpers (shown)
Overall Size: 51.625"L 30.625"D 39"H



Mobile Transport Cage with two solid hinged lids with rod grills welded in the front, sides, and back (shown)
Overall Size: 52.50"L 31.50"D 39"H

These Stainless Steel Mobile Transport Cages with Removable Collection Pans are available with hinged lids or welded tops, welded sides, welded ends or hinged gate on one end. The top, front, sides, and back are available with rod grills, solid panels, or any combination.

Hinged lids open 270 degrees.

The frame is fabricated from die-formed angles, and all joints are heli-arc welded for strength.

The frame and all rod grills have a glass beaded satin finish, removing any discoloration.

The removable glass beaded door has a rod-frame with vertical and horizontal rods and a two-point self-locking latch.

Each of the three (3) removable channel slat floor sections sit on top of the deflector that directs all liquids into the removable collection pans. The cage has two tapered removable collection pans.

The cages are equipped with chrome plated swivel casters with 5" polyolefin wheels, two with brakes. Swivel and axle bearings have seals and zerkl grease fittings.

OPTIONS

12" vertical bumpers
Stainless Steel casters

SSCI reserves the right to change any dimensions.