

# Simplette Basic



Simplify Your Equipment Needs

Simplify your equipment needs with the LABREPCO Simplette Basic, Pipette Controller. The Simplette Basic provides the researcher the highest standards in accuracy and comfort at a lower cost enabling you to focus on what is important, *Advancing Science*.

Compress the suction bellows before attaching the pipet.

Expel residual liquid when using blow-out pipets.

Replace the filter and clean the adapter if necessary after unscrewing the adapter support



Filling and delivery are easily controlled with the small sensitive lever.

**Features:**

- Lightweight ergonomic design.
- Comfortable and simple to use.
- Precise pipetting control.
- Fits entire range of volumetric and graduated pipets, both glass and plastic, from 0.1 to 100 ml.
- Uses replaceable filter.
- Easy to maintain and clean

Cat. No	Color
LAB-100827	blue
LAB-500827	pink



\*\*Specifications are subject to change without prior notice

*\*Pipette tips and racks sold separately.*

*Other products available include:*



Simplette  
Single Channel Pipette



Simplette Multi  
Multi Channel Pipette



Simplette Plus  
Motorized Pipette Filler



**LABREPCO, Inc**  
101 Witmer Road, Suite 700  
Horsham, PA 19044  
USA

Tel: 1-800-521-0754  
Fax: 215-442-9202  
www.labrepco.com  
sales@labrepco.com

**LABREPCO**  
Many Companies • One Solution