

# Simplette Multi



Simplify Your Equipment Needs

Simply turn the plunger button for volume selection

Finger support keeps the instrument in place with minimum user effort

Tip ejector allows convenient one-handed operation

Clear and easy to read display

Dispensing head rotates for optimal pipetting convenience

Durable tip cone materials provide excellent chemical resistance

Universal tip fit provides flexibility in obtaining replacements.



Simplify your equipment needs with the LABREPCO Simplette Multi, Multi-Channel Pipette. The Simplette Multi provides the researcher the highest standards in accuracy and comfort at a lower cost enabling you to focus on what is important, *Advancing Science.*

**Features:**

- Lightweight ergonomic design.
- Easy to calibrate and maintain with calibration tool.
- Ejector works with dispensing head in any position.
- 8 and 12 channel pipettors available for standard 96 well.
- Individual piston/tip cone assemblies allow easy repair and maintenance.
- Compound material tip cone design allows visual seal verification.
- Comfort grip reduces physical stress and discomfort.
- Innovative and robust construction.
- Universal tip fit

Cat. No	Volume Range	Increment	Test Volume	Accuracy (+/-)	Precision (+/-)
<b>8 Channel Pipettes</b>					
LAB-012077	0.5-10 µl	0.1 µl	10 µl 5 µl 1 µl	1.50% 2.50% 4.00%	1.50% 2.50% 4.00%
LAB-022077	5-50 µl	0.5 µl	50 µl 25 µl 5 µl	1.00% 1.50% 3.00%	0.50% 1.00% 2.00%
LAB-042077	50-300 µl	5 µl	300 µl 150 µl 50 µl	0.70% 1.00% 1.50%	0.25% 0.50% 0.80%
<b>12 Channel Pipettes</b>					
LAB-013077	0.5-10 µl	0.1 µl	10 µl 5 µl 1 µl	1.50% 2.50% 4.00%	1.50% 2.50% 4.00%
LAB-023077	5-50 µl	0.5 µl	50 µl 25 µl 5 µl	1.00% 1.50% 3.00%	0.50% 1.00% 2.00%
LAB-043077	50-300 µl	5 µl	300 µl 150 µl 50 µl	0.70% 1.00% 1.50%	0.25% 0.50% 0.80%

Other products available include:

\*\*Specifications are subject to change without prior notice

\*Pipette tips and racks sold separately.



Simplette  
Single Channel Pipette



Simplette Plus  
Motorized Pipette Filler



Simplette Basic  
Pipette Controller



**LABREPCO, Inc**  
101 Witmer Road, Suite 700  
Horsham, PA 19044  
USA

Tel: 1-800-521-0754  
Fax: 215-442-9202  
www.labrepcoco.com  
sales@labrepcoco.com

**LABREPCO**  
Many Companies • One Solution