

Simplette Plus



Accuracy

LABREPCO

Comfort

Simplicity

Simplify Your Equipment Needs

Simplify your equipment needs with the LABREPCO Simplette Plus, Motorized Pipette Filler. The Simplette Plus provides the researcher the highest standards in accuracy and comfort at a lower cost enabling you to focus on what is important, *Advancing Science.*



Separate aspirate and dispense speed controls.

Fully autoclavable nose cone housing simplifies disinfection and decontamination.



Large LCD Display provides visual confirmation of remaining battery charge and speed settings.

Select the desired speed by turning the speed adjustment wheel.

Features:

- Lightweight ergonomic design.
- Comfortable and simple to use.
- Full volume range of 0.1-100 ml.
- Self-regulating aspirator and dispensing speed using circuit feedback.
- Lithium-ion battery offers long usage per charge.
- Low battery indicator.
- Powerful motor enabling quick filling of up to 100ml.
- Replaceable Filter.

Other products available include:



Simplette Single Channel Pipette



Simplette Multi Multi Channel Pipette



Simplette Pipette Controller

Technical Specifications	LAB-239017
Aspirate Speeds	8
Dispense Speeds	8 + Gravity Dispense
Battery	Lithium-Ion
Battery Service Life	More than 15 hours of continuous use
Charging Time	Approx. 2 Hours
Pipette Types	Glass or plastic (0.1 - 100ml), pasteur pipettes
Filter	Hydrophobic

***Specifications are subject to change without prior notice*

LABREPCO, Inc
 101 Witmer Road, Suite 700
 Horsham, PA 19044
 USA

Tel: 1-800-521-0754
 Fax: 215-442-9202
 www.labrepco.com
 sales@labrepco.com

LABREPCO
 Many Companies • One Solution