

Simplette



Simplify Your Equipment Needs

Simply turn the plunger button for volume selection.

Finger support keeps the instrument in place with minimum user effort.

Tip ejector allows convenient one-handed operation.

Ejector collar and tip cone can be removed for easy cleaning and maintenance.

Durable tip cone materials provide excellent chemical resistance.

Universal tip fit provides flexibility in obtaining replacements.

Features:

- Lightweight ergonomic design
- Display clearly reads volume setting
- Volume range from 0.1 µl to 5,000 µl
- Comfort grip reduces physical stress and discomfort to avoid injury from repetitive strain
- Innovative and robust construction
- Universal tip fit

**Pipette tips and racks sold separately.*



Other products available include:



Simplify your equipment needs with the LABREPCO Simplette Single Channel Pipette.

The Simplette Multi provides the operator the highest standards in accuracy and comfort at a lower cost enabling you to focus on what is important,
Advancing Science.

Cat. No	Range	Increment	Volume	Accuracy (+/-)	Precision (+/-)
LAB-500077	0.1-2.5 µl	0.05 µl	2.5 µl 1.25 µl 0.25 µl	2.50% 3.00% 12.00%	2.00% 3.00% 6.00%
LAB-000077	0.5-10 µl	0.1 µl	10 µl 5 µl 1 µl	1.00% 1.50% 2.50%	0.80% 1.50% 1.50%
LAB-080077	2-20 µl	0.5 µl	20 µl 10 µl 2 µl	0.90% 1.20% 3.00%	0.40% 1.00% 2.00%
LAB-020077	5-50 µl	05. µl	50 µl 25 µl 5 µl	0.60% 0.90% 2.00%	0.30% 0.60% 2.00%
LAB-050077	10-100 µl	1 µl	100 µl 50 µl 10 µl	0.80% 1.00% 3.00%	0.15% 0.40% 1.50%
LAB-070077	20-200 µl	1 µl	200 µl 100 µl 20 µl	0.60% 0.80% 3.00%	0.15% 0.30% 1.00%
LAB-060077	100-1000 µl	5 µl	1000 µl 500 µl 100 µl	0.60% 0.70% 2.00%	0.20% 0.25% 0.70%
LAB-011077	1000-5000 µl	50 µl	5000 µl 2500 µl 1000 µl	0.50% 0.60% 0.70%	0.15% 0.30% 0.30%

**Specifications are subject to change without prior notice

LABREPCO, Inc
101 Witmer Road, Suite 700
Horsham, PA 19044
USA

Tel: 1-800-521-0754
Fax: 215-442-9202
www.labrepco.com
sales@labrepco.com

LABREPCO
Many Companies • One Solution