


Security Lock Devices for Upright Freezer Racks

(Patent pending)

Protect your samples from unauthorized access using these security lock devices.



Lock Device	Lock Dimension		Applicable to Rack	Rack Height
	L	H		
DP-0507	4 ³ / ₄ inch 120 mm	6 ³ / ₄ -7 ¹ / ₈ inch 171-181 mm	UFM-1305	6 ² / ₃ inch 169 mm
DP-0509	4 ³ / ₄ inch 120 mm	8 ⁷ / ₈ -9 ⁵ / ₁₆ inch 226-236 mm	UFM-1405	8 ¹³ / ₁₆ inch 224 mm
DP-0511	4 ³ / ₄ inch 120 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UFM-1505	11 inch 280mm
DP-0607	5 ¹ / ₂ inch 140 mm	7 -7 ² / ₅ inch 178-188 mm	UFM-1315	6 ¹⁵ / ₁₆ inch 176mm
DP-0609	5 ¹ / ₂ inch 140 mm	9 ⁵ / ₁₆ -9 ³ / ₄ inch 237-247 mm	UFM-1415	9 ¹ / ₄ inch 235mm
DP-0612	5 ¹ / ₂ inch 140 mm	11 ³ / ₅ -12 inch 295-305 mm	UFM-1515	11 ⁹ / ₁₆ inch 293mm
DP-1107	10 ¹ / ₂ inch 259 mm	6 ¹¹ / ₁₆ -7 ¹ / ₈ inch 171-181 mm	UFM-2305	6 ² / ₃ inch 169 mm
DP-1109	10 ¹ / ₂ inch 259 mm	8 ⁷ / ₈ -9 ⁵ / ₁₆ inch 226-236 mm	UFM-2405	8 ¹³ / ₁₆ inch 224 mm
DP-1111	10 ¹ / ₂ inch 259 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UFM-2505	11 inch 280mm
DP-1210	11 ² / ₃ inch 290 mm	9 ¹⁵ / ₁₆ -10 ² / ₅ inch 253-263 mm	UFLB-22	10 ³ / ₁₆ inch 259mm
DP-1307	11 ¹ / ₄ inch 298 mm	7 -7 ² / ₅ inch 178-188 mm	UFHT-23	6 ¹⁵ / ₁₆ inch 176mm
			UFM-2315	6 ¹⁵ / ₁₆ inch 176mm
DP-1309	11 ³ / ₄ inch 298 mm	9 ⁵ / ₁₆ -9 ³ / ₄ inch 237-247 mm	UFHT-24	9 ¹ / ₄ inch 235mm
			UFM-2415	9 ¹ / ₄ inch 235mm
DP-1312	11 ¹ / ₄ inch 298 mm	11 ³ / ₅ -12 inch 295-305 mm	UFHT-25	11 ⁹ / ₁₆ inch 293mm
			UFM-2515	11 ⁹ / ₁₆ inch 293mm
DP-1607	15 ¹¹ / ₁₆ inch 395 mm	6 ¹¹ / ₁₆ -7 ¹ / ₈ inch 170-180 mm	UF-332	6 ⁵ / ₈ inch 168 mm
			UF-323	6 ⁷ / ₈ inch 168 mm
DP-1609	15 ⁹ / ₁₆ inch 395 mm	8 ¹⁵ / ₁₆ -9 ⁵ / ₁₆ inch 227-237 mm	UF-342	8 ⁷ / ₈ inch 225 mm
DP-1610	15 ⁹ / ₁₆ inch 395 mm	9 ² / ₅ -9 ⁹ / ₁₆ inch 239-249 mm	UF-333	9 ³ / ₈ inch 238 mm
DP-1611	15 ⁹ / ₁₆ inch 395 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UF-352	11 inch 280mm
DP-1613	15 ⁹ / ₁₆ inch 395 mm	12 ³ / ₄ -13 ¹ / ₈ inch 323-333 mm	UF-362	12 ² / ₈ inch 321mm
DP-1707	15 ¹¹ / ₁₆ inch 398 mm	6 ³ / ₄ -7 ¹ / ₈ inch 171-181 mm	UFM-3305	6 ² / ₃ inch 169 mm
DP-1709	15 ¹¹ / ₁₆ inch 398 mm	8 ⁷ / ₈ -9 ⁵ / ₁₆ inch 226-236 mm	UFM-3405	8 ¹³ / ₁₆ inch 224 mm
DP-1711	15 ¹¹ / ₁₆ inch 398 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UFM-3505	11 inch 280mm

(Continued on next page)

(Continued from last page)

Security Lock Devices for Upright Freezer Racks

Lock Device	Lock Dimension		Applicable to Rack	Rack Height
	L	H		
DP-1809	17 ⁷ / ₁₆ inch 446 mm	8 ⁷ / ₁₆ -9 ⁵ / ₁₆ inch 226-236 mm	UFTB-54	8 ¹³ / ₁₆ inch 224 mm
		9 ¹ / ₄ -10 ¹ / ₁₆ inch 235-255 mm	UFDP-55	9 ³ / ₁₆ inch 233 mm
			UFMP-509	9 ⁹ / ₁₆ inch 233 mm
DP-1810	17 ² / ₅ inch 443 mm	9 ¹⁵ / ₁₆ -10 ⁵ / ₁₆ inch 252-262 mm	UFLB-32	10 ³ / ₁₆ inch 259mm
DP-1907	17 ¹⁵ / ₁₆ inch 456 mm	7 -7 ² / ₈ inch 178-188mm	UFHT-33	6 ¹⁵ / ₁₆ inch 176mm
			UFM-3315	6 ¹⁵ / ₁₆ inch 176mm
DP-1909	17 ¹⁵ / ₁₆ inch 456 mm	9 ⁵ / ₁₆ -9 ³ / ₄ inch 237-247mm	UFHT-34	9 ¹ / ₄ inch 235mm
			UFM-3415	9 ¹ / ₄ inch 235mm
DP-1911	17 ³ / ₈ inch 446 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292mm	UFDP-56	11 inch 280mm
			UFTB-55	11 inch 280mm
			UFMP-511	11 inch 280mm
DP-1912	17 ¹⁵ / ₁₆ inch 456 mm	11 ³ / ₈ -12 inch 295-305 mm	UFHT-35	11 ⁹ / ₁₆ inch 293mm
			UFM-3515	11 ⁹ / ₁₆ inch 293mm
DP-2107	20 ⁷ / ₈ inch 530 mm	6 ¹¹ / ₁₆ -7 ¹ / ₈ inch 170-180 mm	UF-432	6 ⁵ / ₈ inch 168 mm
			UF-423	6 ⁵ / ₈ inch 168 mm
DP-2109	20 ⁷ / ₈ inch 530 mm	8 ¹⁵ / ₁₆ -9 ⁵ / ₁₆ inch 227-237 mm	UF-442	8 ⁷ / ₈ inch 225 mm
DP-2110	20 ⁷ / ₈ inch 530 mm	9 ² / ₅ -9 ⁴ / ₅ inch 239-249 mm	UF-433	9 ³ / ₈ inch 238 mm
DP-2111	20 ⁷ / ₈ inch 530 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UF-452	11 inch 280mm
DP-2112	21 ¹ / ₈ inch 542 mm	12 ³ / ₄ -13 ¹ / ₈ inch 323-333 mm	UF-443	12 ³ / ₈ inch 321 mm
DP-2113	20 ⁷ / ₈ inch 530 mm	12 ³ / ₄ -13 ¹ / ₈ inch 323-333 mm	UF-462	12 ⁵ / ₈ inch 321 mm
DP-2207	21 ¹ / ₈ inch 536 mm	6 ¹ / ₄ -7 ¹ / ₈ inch 171-181 mm	UFM-4305	6 ² / ₅ inch 169 mm
DP-2209	21 ¹ / ₅ inch 538 mm	9 ¹ / ₄ -10 ¹ / ₁₆ inch 235-255 mm	UFDP-65	9 ⁹ / ₁₆ inch 233 mm
			UFMP-609	9 ⁹ / ₁₆ inch 233 mm
		8 ⁷ / ₈ -9 ⁵ / ₁₆ inch 226-236 mm	UFTB-64	8 ¹³ / ₁₆ inch 224 mm
			UFM-4405	8 ¹³ / ₁₆ inch 224 mm
DP-2211	21 ¹ / ₅ inch 538 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UFDP-66	11 inch 280mm
			UFTB-65	11 inch 280mm
			UFMP-611	11 inch 280mm
			UFM-4505	11 inch 280mm
DP-2410	23 ³ / ₅ inch 600 mm	9 ¹⁵ / ₁₆ -10 ² / ₅ inch 253-263 mm	UFLB-42	10 ³ / ₁₆ inch 259mm
DP-2609	24 ⁴ / ₅ inch 631 mm	9 ¹ / ₄ -9 ³ / ₄ inch 235-245 mm	UFMP-709	9 ³ / ₁₆ inch 233 mm
DP-2611	24 ⁴ / ₅ inch 631 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UFMP-711	11 inch 280mm
DP-2707	25 ¹⁵ / ₁₆ inch 659 mm	6 ¹ / ₁₆ -7 ¹ / ₈ inch 170-180 mm	UF-532	6 ⁵ / ₈ inch 168mm
DP-2709	25 ¹⁵ / ₁₆ inch 659 mm	8 ¹⁵ / ₁₆ -9 ⁵ / ₁₆ inch 227-237 mm	UF-542	8 ⁷ / ₈ inch 225mm
DP-2710	25 ¹⁵ / ₁₆ inch 659 mm	9 ⁵ / ₁₆ -9 ¹⁵ / ₁₆ inch 243-253 mm	UF-533	9 ¹ / ₂ inch 241mm
DP-2711	25 ¹⁵ / ₁₆ inch 659 mm	11 ¹ / ₈ -11 ¹ / ₂ inch 282-292 mm	UF-552	11 inch 280mm
DP-2713	25 ¹⁵ / ₁₆ inch 659 mm	12 ³ / ₄ -13 ¹ / ₈ inch 323-333 mm	UF-543	12 ⁵ / ₈ inch 321mm

- *No tools are required for installing lock devices.
- *Locking devices are made of durable Stainless Steel.
- *Protects against unwanted access to freezer samples.
- *Select lock part number for your specific Upright Freezer Rack.