

Lab Animal Research (LAR) Carts — 52.05

LAR Feed Carts are designed to store and transport animal feed. Stainless steel finish provides superior corrosion resistance, while compact footprint makes it easy to maneuver in and around tight lab spaces. Clear totes provide easy identification of tote contents.

Lab Security Cart provides secure storage for animal cage, lab instrument, or hazardous chemical transport. Solid stainless steel shelf with raised ship's edge can be used as a work surface and can contain unwanted spills.

- Super Erecta Shelf® construction makes units easy to assemble with absolute rigidity, and makes relocation of shelves quick and simple.
- Wire Shelves: Bright, modern, sanitary appearance. Open construction minimizes dust accumulation, maximizes visibility and air circulation.

Width (in.)	Length (in.)	Height (in.)	Rating ¹ (lbs.)	Description	Approx. Pkd. Wt. (lbs.)	Cat. No.	List Price Each					
18	457	24	610	38	965	375	171	Light-duty feed cart One tote capacity/shelf 4" (102mm) resilient casters ²	42	18.7	LAR1	669.50
24	610	36	914	39	990	600	273	Heavy-duty feed cart Two tote capacity/shelf 5" (127mm) polyurethane casters ²	65	29.2	LAR2	1,037.00
24	610	30	760	39	990	600	273	Stainless steel security cart with flip-up door 5" (127mm) polyurethane casters	85	38.6	SECMLAB³	1,246.50

¹Cart load rating listed is based on an evenly distributed load.

²Cat. No. LAR1 features four L4D series swivel casters with resilient rubber tread. Cat. No. LAR2 features four 5MP series swivel casters with polyurethane tread.

³Lab Security Cart features two stainless steel wire shelves, one stainless steel security module, and one stainless steel solid shelf. Cart comes standard with two brake casters (Cat. No. 5MPB) and two swivel casters (Cat. No. 5MP).



Feed Cart
(One Tote Capacity)



Feed Cart
(Two Tote Capacity)



Lab Security Cart
(shown with optional accessories)



Clear Tote Boxes (shown with Snap-on Cover and Cardholder Accessories)

Refer to page 183 for information on Metro's Polymer LAR Feed Carts.

Accessories

Durable, clear totes provide clean and secure storage for animal feed. Totes are injection-molded from clear, polypropylene material, are stackable, and can be combined with clear snap-on covers and cardholders for security and identification of contents.

Clear Tote Box

(in.)	Top		Outside Dimensions Length/Width/Height		Bottom		Inside Dimensions Length/Width/Height		Bottom	(mm)	Cat. No.	Clear Polypropylene	List Price Each
	(in.)	(mm)	(in.)	(mm)	(in.)	(mm)	(in.)	(mm)					
22.5x17.5x12	572	445	305	21x16	533x406	20.2x15.1x11.5	513x384x292	TB93120CLR	55.00				

Note 1: All outside dimensions have an overall tolerance of (+/-) .050". All inside dimensions have an overall tolerance of (+/-) .100"

Note 2: Tote not recommended for freezer environments (temperatures lower than -20°F [-29°C]).

Clear Snap-on Cover

Length (in.)	(mm)	Outside Dimensions		Height (in.)	(mm)	Cat. No.	List Price Each
		Width (in.)	(mm)				
22.5	572	17.5	445	1.25	32	CO9300CLR	16.00

Tote Box Cardholders

Material	Inside Dimensions				Cat. No.	List Price Each
	Length (in.)	(mm)	Height (in.)	(mm)		
Vinyl	4.5	114.3	8	203	OP2501CLR	3.50
Vinyl	2 ¹¹ / ₁₆	60.6	5	127	OP2535CLR	3.50

Note: Pricing based on quantity ordered. Please consult your Metro representative.

Cover Holder

Stainless steel wire cover holder provides a convenient storage option for tote box covers when not in use. Accessory easily mounts on any Super Erecta® based wire shelf.

Width (in.)	(mm)	Length (in.)	(mm)	Depth (in.)	(mm)	Approx. Carton Wt. (lbs.)	(kg)	Cat. No.	List Price Each
18 ¹ / ₂	470	18 ¹ / ₄	464	4	101	6	2.7	LAR-CH	69.00



Cover Holder (shown mounted on Cart with one Tote Cover)