



GLACIER Freezers -86°C

Quality, Dependability and Priced Right.



Multi-point gasket seal

Vacuum Panel Insulation
High density foam produces greater thermal insulation properties (NU-9728 only)

Vacuum Relief Valve

Noise Reduction
Condenser fan design and noise reduction insulation suppresses sound as low as 47 dba.

Remote alarm
contact and optional communication port

Universal key door lock

Access Ports
Multiple access ports allow for the use of back-up systems, temperature recording, and monitoring systems.

Efficient Refrigeration System
Twin Compressors operate high and low stage HCFC-free refrigerants. Heat generated is dispersed in multiple stages of cooling. This will prolong the life of the compressor.

Easy access
washable pre-filter

Optional circular chart
recorder mounting space

Inner Door
Insulated and gasket door with claw latch

Claw Door Latch
Durable, Easy to use Claw Door Latch with integrated pad lock

Ergonomic Eye-Level Controls

Insulated Inner Doors minimize temperature loss during routine door openings and reduce frost build-up.

Glacier NU-9483GA or GC	Int. Volume 17.1 ft ³	Electrical GA : 115V/60Hz GC : 220V/60Hz	Configuration Racks - 4 columns / 4 Rows # of 2+ 384
Glacier NU-9668	Int. Volume 23.5 ft ³	Electrical GC : 220V/60Hz	Configuration Racks - 6 columns / 4 Rows # of 2+ 480
Glacier NU-9728	Int. Volume 25.7 ft ³	Electrical GA : 115V/60Hz	Configuration Racks - 6 columns / 4 Rows # of 2+ 576

For over 40 years, NuAire has provided the life science and medical research communities with valued products that are highly dependable, effective and efficient.

Call or Email NuAire's Local Sales Representative

- LabrepcO, Inc. - 800-521-0754 - sales@labrepcO.com -