

basic

fixed laboratory system



CORE SUPPORT

- FULL SEFA 8 Compliance
- Removable back panels for easy service access
- Epoxy, Phenolic, Stainless work top options
- Full array of electrical and plumbing service solutions.

forte™

adaptable laboratory system



The Forte™ system's tables provide simplicity and strength at it's best.

- 1500 lb. weight load capacity
- Height adjustable worktops: 30"-36" in 1" increments
- Single or double sided configurations
- Integrated plug and play services

ergolab™

mobile laboratory system



Mobile or Stationary...
...the best of both worlds

- 2600 lb. bench load capacity
- Vertically adaptable from 36" to 77"
- Caster-based design with full leveling capability
- Comprehensive array of ergonomic accessories

bravo

laboratory benching system



Winner: Lab of the Year 2012

- 1500 lb. weight load capacity
- Vertically adaptable from 36" to 84"
- Integrated plug and play services
- Comprehensive array of ergonomic accessories