



NuAire Centrifuges

NuAire is the new optimal choice in Centrifugation. The NuAire centrifuge product line now extends from Mini-centrifuges, Micro-centrifuges, and General Purpose Multi-Application Centrifuges. Each centrifuge is designed for simplicity and functionality to everyday use.

Join the NuAire family with a NuAire Centrifuge today.



Crystal

Mini-centrifuges


Capacity	16 x 0.2 ml (Rectangular) 6 x 1.5 ml (Round)
Max Speed	
Pink	4000 rpm (900 xg)
White	6,000 rpm (2,000 xg)
Blue	10,000 rpm (5,000 xg)
Dims (W x D x H)	155 x 140 x 115 mm 6 x 5.5 x 4.5 in.



Breeze

Microcentrifuges

Capacity	24 x 1.5 ml 16 x 0.2 ml
Max Speed	
Ventilated	15,000 rpm (21,200 xg)
Refrigerated	13,500 rpm (17,135 xg)
Dims (W x D x H)	
Ventilated	240 x 350 x 190 mm 9.45 x 13.8 x 7.48 in.
Refrigerated	277 x 450 x 247.5 mm 10.9 x 17.75 x 9.75 in.



NuWind

Multi-Application Centrifuges

Capacity	92 x 5 ml 76 x 10 ml 44 x 15 ml Conical 20 x 50 ml Conical 16 Microplate
Max Speed	
Swing-Out	5,100 rpm (5,350 xg)
Angled	18,000 rpm (30,980 xg)
Dims (W x D x H)	
Ventilated	457 x 604 x 347 mm 18.1 x 23.9 x 13.7 in.
Refrigerated	675 x 630 x 322 mm 26.7 x 25 x 12.7 in.



Your Local NuAire Representative

NuAire, Inc. | 2100 Fernbrook Lane | Plymouth, MN 55447 | U.S.A. | 763.553.1270 | www.nuaire.com

