



# COMPETITION

## Ventilated Model



**NuAire**  
**NU-C200V**

**VS**



**Eppendorf**  
**5810**

Best Products. | Best Performance. | Best Protection.



# COMPETITION

## Ventilated Model





Best Products. | Best Performance. | Best Protection.



# COMPETITION

## Ventilated Model




Best Products. | Best Performance. | Best Protection.



# COMPETITION

## Ventilated Model

Size	NuAire NU-C200V	Eppendorf 5810
Dimensions (HxWxD)	13.7x18.1x22	13x21x21
Front view		
Top view		

Best Products. | Best Performance. | Best Protection.



# COMPETITION

## Ventilated Model

Specs		NuAire NU-C200V		Eppendorf 5810	
Max volume		4x500ml		4x500ml	
Max speed	S-out	4,500 RPM / 4,188 xg		4,000 RPM / 3,250 xg	
	Angle	15,300 RPM / 22,380 xg		14,000 RPM / 20,800 xg	
Max volume	5ml	92 tubes	8% less	100 tubes	
	10ml	76 tubes	18% more	64 tubes	
	15ml con	44 tubes	8% less	48 tubes	
	50ml con	20 tubes	same	20 tubes	
	Plates	16 standard	same	16 standard	

Best Products. | Best Performance. | Best Protection.



# COMPETITION

## Ventilated Model



- NuAire C200V:

- + • Modern design
- Easy / Quick to service
- Efficient ventilation, keeps samples at ambient temperature
- Touch screen
- Higher speeds
- No tools to change rotors
- • Reputation in centrifugation

- Eppendorf 5810:

- + • More market shares
- Reliable and great quality
- Large selection of accessories
- Brand reputation
- • Very expensive
- Old design
- Need tool to change rotor



# COMPETITION

## Ventilated Model



**NuAire**  
NU-C200V



**Eppendorf**  
5810

### UNIQUE FEATURES:

- Control Center to provide constant monitoring, easy servicing, no down time
- Diagnostic via software
- Preventive maintenance alerts
- Max speed of 4,500RPM (S-out) / 15,300RPM (Angle)
- Motor can be removed from the top for servicing purpose
- Excellent ventilation with no recirculation of air
- Touch screen display
- Interchangeability of rotors with no tool
- Power Factor Control

Best Products. | Best Performance. | Best Protection.

NuAire, Inc. | 2100 Fernbrook Lane | Plymouth, MN 55447 | U.S.A | ph: 763.553.1270 | fx: 763.553.0459 | tf: 800.328.3352 | [www.nuaire.com](http://www.nuaire.com)





# COMPETITION

## Refrigerated Model



**NuAire**  
**NU-C200R**

**VS**



**Eppendorf**  
**5810R**

Best Products. | Best Performance. | Best Protection.





# COMPETITION

## Refrigerated Model





Best Products. | Best Performance. | Best Protection.



# COMPETITION

## Refrigerated Model





Best Products. | Best Performance. | Best Protection.



# COMPETITION

## Refrigerated Model

Size	NuAire Nu-C200R	Eppendorf 5810R
Dimensions (HxWxD)	12.7x26.7x23	13x28x21
Front view		
Top view		

Best Products. | Best Performance. | Best Protection.



# COMPETITION

## Refrigerated Model

Specs		NuAire NU-C200R		Eppendorf 5810R	
Max volume		4x500ml		4x500ml	
Max speed	S-out	5,100 RPM / 5,350 xg		4,000 RPM / 3,250 xg	
	Angle	18,000 RPM / 30,980 xg		14,000 RPM / 20,800 xg	
Max volume	5ml	92 tubes	8% less	76 tubes	
	10ml	76 tubes	18% more	56 tubes	
	15ml con	44 tubes	8% less	36 tubes	
	50ml con	20 tubes	same	16 tubes	
	Plates	16 standard	same	6 standard (separate rotor)	

Best Products. | Best Performance. | Best Protection.



# COMPETITION

## Refrigerated Model



- NuAire C200R:

- + • Modern design
- Easy to service
- Touch screen
- Much higher speeds
- Constant monitoring
- Complete package w/ NuAire items
- Extremely silent
- • Reputation in centrifugation
- Max volume of 4x500ml

- Eppendorf 5810R:

- + • More market shares
- Reliable and great quality
- Large selection of accessories
- Brand reputation
- Very good refrigeration
- 4x750ml in the same shell
- • Very expensive
- Old design
- Need tool to change rotor

Best Products. | Best Performance. | Best Protection.



# COMPETITION

## Refrigerated Model



**NuAire**  
NU-C200R



**Eppendorf**  
5810R

### UNIQUE FEATURES:

- Control Center to provide constant monitoring, easy servicing, no down time
- Diagnostic via software
- Preventive maintenance alerts
- Max speed of 5,100RPM (S-out) / 18,000RPM (Angle)
- Motor can be removed from the top for servicing purpose
- Excellent ventilation with no recirculation of air
- Touch screen display
- Interchangeability of rotors with no tool
- Power Factor Control

Best Products. | Best Performance. | Best Protection.

NuAire, Inc. | 2100 Fernbrook Lane | Plymouth, MN 55447 | U.S.A | ph: 763.553.1270 | fx: 763.553.0459 | tf: 800.328.3352 | [www.nuaire.com](http://www.nuaire.com)