



Introducing the most energy efficient ultra-low temperature freezer on the market.¹

ENERGY STAR® Certified and independently tested by a nationally recognized testing laboratory.

- 7.30 kWh/day, Steady State Energy Consumption
- 7.87 kWh/day, Daily Energy Consumption

Model MDF-DU702VH-PA



VIP® ECO

Your Life's Work, Our Brand.

The Panasonic Healthcare -86°C VIP ECO ultra-low temperature freezer delivers the industry's most energy efficient performance. Because our unique variable differential cascade refrigeration system uses all natural refrigerants, you can achieve the quality of cold you need with the global sustainability you desire.

Fast recovery. Impeccable uniformity. High ambient tolerance. Learn more about the Panasonic Healthcare VIP ECO 25.7 cu.ft. (729 L) upright ultra-low temperature freezer. Designed, tested and field-proven for stable, reliable storage of biologicals.

www.panasonic-healthcare.com/us/biomedical/vip-eco



¹Based on independent third-party testing at time of publication. ENERGY STAR test results for submitted products can be compared for performance across the competitive market.

Results are published on the ENERGY STAR website www.energystar.gov. See Panasonic Certification Number 4787624501.