

# Top 3 Reasons to Store Your Vaccines in a Clinical Unit:

1. The CDC recommends storing all vaccines (or other temperature sensitive medications) in stand-alone, pharmacy grade/ purpose-built units.



2. All refrigerators come standard with microprocessor temperature controllers for excellent temperature accuracy, digital temperature displays, audible & visual hi/low temperature alarms, and internal fans for even temperature distribution throughout. (Refrigerators Only)



3. All units maintain a strict temperature for the storage of refrigerated (+35°F to 46°F) and frozen (-58F to +5°F) vaccines.

[Click Here for Clinical Refrigerators](#)

[Click Here for Clinical Freezers](#)

