

Our guide to help you find the right size refrigerator or freezer for your vaccine storage

Vaccine Storage requiring a refrigerator

Click on model #
for more information

Max Doses	Minimum Cubic Feet	Recommended Refrigerators
0-150	3	LHP-3-URG-PH (glass door) LHP-3-UR-PH (solid door)
151-399	5	LHP-5-URBSS-PH (solid door) LHP-5-URBG-PH (glass door) LHP-5-URBGPH-ADA (glass door) LHP-5-URB-PH (solid door)
400-700	11 to 12	LHP-11-HG-PH (glass door) LHU-12-HG-PH (glass door)
701-900	23	LHP-23-SDSS-PH (solid door) LHU-23-HG-PH (glass door)
901-2,000	49	LHP-49-SDSS-PH (solid doors) LHU-49-HG-PH (glass doors)



Vaccine Storage Requiring a Freezer

Max Doses	Minimum Cubic Feet	Recommended Refrigerators
0-51	2	LHP-2-FM-PH
51-200	5	LHP-5-UFMBSS-PH LHP-5-UFMB-PH LHP-5-UFMBPH-ADA
201-600	20	LHE-20-FM

Can't find the size you need?
Give us a call at 1-800-521-0754



CliniCool
LABRepCo



Ask about our Hydrocarbon
"HC" Refrigerant

