



LAAB Elite – Active Package

\$7,200

Take your lab anywhere with the Lab Animal Anesthesia Bench (LAAB) – Elite. The Link7 is teamed up with a portable oxygen source, a magnetic workspace and two shelves for anesthesia accessory storage. Whether your anesthesia system is used by multiple departments or simply needs to be moved from room to room, this system will meet your needs. System comes complete including the cart with magnetic workspace, a dual e-cylinder manifold with O2 hose, the Link 7 system, vaporizer, Posi-vac nosecone, and hooded induction chamber. WAG filters are also included to adsorb the waste gas that will be actively managed through the Link 7 system.

Package includes:

- 78930337 MOBILE LAAB ELITE SYSTEM (includes Link 7, EX3 Vaporizer-6psi, Dual e-cyl manifold, cart)
- 78930339 MOBILE LAAB POLE MOUNT OPTION
- 78040735 ANTI SPILL ISOFLURANE/FLORANE
- 78914728 POSI-VAC NRB SYSTEM - MOUSE/SMALL RAT
- 78933385 HOODED INDUCTION CHAMBER PVC 12X6X8
- 78909583 WAG CANISTER HOLDER
- 78909457 WAG ACTIVATED CHARCOAL FILTER CA4