

Job _____



D. For Autoclave (Mobile) Applications (cont'd.)

Autoclavable Cage Racks

- Preconfigured, autoclavable cage racks available in two footprints — 18x60 (457x1524mm) and 24x60 (610x1524mm).
- Units are 69" (1753mm) high, including casters. **Casters must be ordered separately.**
- Carts come standard with swedged posts and aluminum split sleeves.

Cat. No.	Width (in.)	Length/Height (mm)	Approx. Pkd. Wt. (lbs.) (kg)		Description
LC5S1860	18x60x69	457x1524x1753	136	61.7	5-shelf, 18" (457mm) wide cage rack
LC5S2460	24x60x69	610x1524x1753	166	75.3	5-shelf, 24" (610mm) wide cage rack

Note: To build a cage rack according to specific requirements, follow instructions (D) on page 7. Refer to Autoclave guidelines, below, for important checklist to follow

Autoclave Cart Covers

Cat. No.	Width (in.)	Length/Height (mm)	Approx. Pkd. Wt. (lbs.) (kg)		Description
LC1860AC	20x62x61.5	508x1575x1562	6.0	2.7	Autoclavable cart cover for 18" (457mm) wide Cage Rack
LC2460AC	26x62x61.5	660x1575x1562	7.0	3.2	Autoclavable cart cover for 24" (610mm) wide Cage Rack

Note: Autoclave cart covers are rated for approximately 70 cycles per cover.

Autoclavable Lab Carts

- Preconfigured, autoclavable Lab carts are available in two footprints — 18x24 (457x610mm) and 24x36 (610x914mm).
- Chose between 2- and 3-shelf models
- Units are 39" (991mm) high, including casters. **Caster must be ordered separately.**
- Carts come standard with swedged posts and aluminum split sleeves.

Cat. No.	Width (in.)	Length/Height (mm)	Approx. Pkd. Wt. (lbs.) (kg)		Description
LC2S1824	18x24x39	457x610x991	25.0	11.3	2-shelf, 18" (457mm) wide Lab cart
LC3S1824	18x24x39	457x610x991	37.5	17.0	3-shelf, 18" (457mm) wide Lab cart
LC2S2436	24x36x39	610x914x991	46.5	21.1	2-shelf, 24" (610mm) wide Lab cart
LC3S2436	24x36x39	610x914x991	65.5	29.7	3-shelf, 24" (610mm) wide Lab cart

Note: To build a lab cart according to specific requirements, follow instructions (D) on page 7. Refer to Autoclave guidelines, below, for important checklist to follow.

Autoclave Cart Covers

Cat. No.	Width (in.)	Length/Height (mm)	Approx. Pkd. Wt. (lbs.) (kg)		Description
LC1824AC	20x26x33.5	508x660x851	3.0	1.4	Autoclavable cart cover for 18" (457mm) wide Cage Rack
LC2436AC	26x38x33.5	660x965x851	3.5	1.6	Autoclavable cart cover for 24" (610mm) wide Cage Rack

Note: Autoclave cart covers are rated for approximately 70 cycles per cover.



LC5S2460
(casters ordered separately)



LC5S2460 with autoclavable cart cover
(casters ordered separately)



LC2S1824
(casters ordered separately)

AUTOCLAVE APPLICATION GUIDELINES

When stainless solid shelving, as part of a mobile shelving unit, is used in an autoclave, the following Metro product guidelines are recommended to ensure optimal product performance and longevity:
 Stainless, swedged posts with aluminum post caps must be utilized [part numbers 33/54/63UPS-SW]
 Aluminum split sleeves [part number 9986S] with stainless rings should be used to mount solid stainless shelves.
 Select either hi-temp phenolic casters or hi-temp nylon casters for application. These casters are designed specifically for autoclave environments.
 Bumpers should not be used as they are not designed to withstand autoclave temperature ranges

All Metro Catalog Sheets are available on our website: www.metro.com



InterMetro Industries Corporation

North Washington Street • Wilkes-Barre, PA 18705 • 570.825.2741
 Fax: 800.638.9263 (East Coast/Canada) • Fax: 800.638.3292 (West Coast)

FOR PRODUCT INFORMATION/CUSTOMER SERVICE:
 U.S./Canada/Latin America: 1.800.992.1776 • Europe: +31.76.587.7550
 Asia/Pacific: +65.6567.8003 • Middle East/Africa: +971.4.811.8286

Information and specifications are subject to change without notice. Please confirm at time of order

L03-259
02/12
Printed in U.S.A.

Copyright © 2012 InterMetro Industries Corp.

