


# Vertical Rack



H:  $\frac{9}{16}$ " -  $1 \frac{3}{4}$ " (14.4-44.32mm)

## For Matrix

**500 $\mu$ l Screw-Top Latch Rack (H=1 $\frac{3}{4}$ inch / 44.32mm)**

**0.75ml Latch Rack (H=1 $\frac{1}{2}$ inch / 39.37mm)**

## For Micronic

**Roborack-96 with Low Cover (H=1 $\frac{5}{8}$ inch / 41.2mm)**

## For Fluidx


**384-Format Rack for Square Tube (H= $\frac{9}{16}$ inch / 14.4mm, each compartment holds 3 racks.)**

**Low Profile 96-Format Rack for all 0.5ml and 0.75ml Tube with Tall Profile Lid**

**(H=1 $\frac{3}{4}$ inch / 43.87mm)**

**96-Format Rack for 1.0ml Tube with Low Profile Lid (H=1 $\frac{3}{4}$ inch / 43.9mm)**

**96-Format Rack for Tube with Protective Jacket with Standard Profile Lid (H=1 $\frac{3}{4}$ inch / 44mm)**



- ◆ Made of corrosion resistant stainless steel.
- ◆ Locking rods are included.

Part Number	Dimension H x D x W	Compartment Size h x d x w	Capacity	Lock Device Option*	Shipping Weight lb(Kg)
CFDP-9	16 $\frac{3}{8}$ x 3 $\frac{5}{8}$ x 5 $\frac{1}{2}$ inch 423 x 92 x 140 mm	1 $\frac{13}{16}$ x 3 $\frac{7}{16}$ x 5 $\frac{3}{16}$ inch 45.5 x 88 x 131.5 mm	9	DR-15	3.1(1.39)
CFDP-10	18 $\frac{7}{16}$ x 3 $\frac{5}{8}$ x 5 $\frac{1}{2}$ inch 469 x 92 x 140 mm	1 $\frac{13}{16}$ x 3 $\frac{7}{16}$ x 5 $\frac{3}{16}$ inch 45.5 x 88 x 131.5 mm	10	DR-17	3.3(1.48)
CFDP-11	20 $\frac{1}{4}$ x 3 $\frac{5}{8}$ x 5 $\frac{1}{2}$ inch 515 x 92 x 140 mm	1 $\frac{13}{16}$ x 3 $\frac{7}{16}$ x 5 $\frac{3}{16}$ inch 45.5 x 88 x 131.5 mm	11	DR-19	3.5(1.61)
CFDP-12	22 $\frac{1}{8}$ x 3 $\frac{5}{8}$ x 5 $\frac{1}{2}$ inch 562 x 92 x 140 mm	1 $\frac{13}{16}$ x 3 $\frac{7}{16}$ x 5 $\frac{3}{16}$ inch 45.5 x 88 x 131.5 mm	12	DR-22	3.8(1.73)
CFDP-13	23 $\frac{15}{16}$ x 3 $\frac{5}{8}$ x 5 $\frac{1}{2}$ inch 608 x 92 x 140 mm	1 $\frac{13}{16}$ x 3 $\frac{7}{16}$ x 5 $\frac{3}{16}$ inch 45.5 x 88 x 131.5 mm	13	DR-23	4.1(1.86)
CFDP-14	25 $\frac{3}{4}$ x 3 $\frac{5}{8}$ x 5 $\frac{1}{2}$ inch 654 x 92 x 140 mm	1 $\frac{13}{16}$ x 3 $\frac{7}{16}$ x 5 $\frac{3}{16}$ inch 45.5 x 88 x 131.5 mm	14	DR-25	4.4(1.99)
CFDP-15	27 $\frac{9}{16}$ x 3 $\frac{5}{8}$ x 5 $\frac{1}{2}$ inch 700 x 92 x 140 mm	1 $\frac{13}{16}$ x 3 $\frac{7}{16}$ x 5 $\frac{3}{16}$ inch 45.5 x 88 x 131.5 mm	15	DR-27	4.7(2.12)
CFDP-16	29 $\frac{7}{16}$ x 3 $\frac{5}{8}$ x 5 $\frac{1}{2}$ inch 747 x 92 x 140 mm	1 $\frac{13}{16}$ x 3 $\frac{7}{16}$ x 5 $\frac{3}{16}$ inch 45.5 x 88 x 131.5 mm	16	DR-29	5.0(2.25)

\* Lock device is ordered separately.