



Silver Series Undercounter Pharmacy/Vaccine Refrigerators

Our CliniCool® Silver Series undercounter pharmacy/vaccine refrigerators are ideal for clinical labs, medical offices and pharmacies with a limited space and budget. These units offer superior temperature performance for safe storage of vaccines while meeting federal and state guidelines.



LHP-3-UR-PH

Standard Product Features:

- Cycle defrost
- Microprocessor temperature controller with battery back-up (batteries included)
- Adjustable operating temperature range: +36°F to +46°F (+2°C to +8°C)
- Freestanding thermometer with 2 year certificate of calibration
- Digital temperature display
- Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable
- Audible and visual high and low temperature alarms, with battery back up
- Remote alarm contacts
- 2 probes; 1 in air 1 in glycol bottle
- Four adjustable shelves (two large, two small) with guard rail on back to prevent product from being pushed against the back wall
- Bottom footed shelf to allow air circulation
- Back wall evaporator cover
- Forced draft circulation
- Keyed door lock
- LED light
- Magnetic door gasket for positive seal
- Leveling legs
- Probe access port (3/8")
- Powder coat exterior
- High density urethane foam cabinet and door insulation
- HFC-free refrigerant (R600a), foam insulation and packaging
- Pharmacy refrigerator/freezer toolkit and temperature logs
- Pyxis®, Omnicell® and AcuDose RX® compatible
- Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage.
- Warranty: 2 years parts and labor

Optional Accessories:

- Solid Window Film *

**For glass door models only. This window film is a thin, solid piece of laminate that covers the glass door which prevents visibility or room light into the cabinet.*



One Company • Many Solutions

Glass Door Models

Solid Door Models

LHP-5-URG-PH

LHP-5-URBSS-PH

LHP-5-URB-PH



LHP-3-URG-PH

LHP-2-URG-PH

LHP-2-UR-PH

LHP-3-UR-PH



Hydrocarbon Refrigerant used in all CliniCool Silver Series undercounter refrigerators

Model	LHP-2-UR-PH	LHP-2-URG-PH	LHP-3-UR-PH	LHP-3-URG-PH	LHP-5-URG-PH	LHP-5-URB-PH	LHP-5-URBG-PH**	LHP-5-URBSS-PH
Weight lbs	40	40	68	68	93	97	97	100
Overall Height in.	21 1/4*	21 1/4*	28 7/8*	28 7/8*	32 1/8*	33 3/8*	33 3/8*	33 3/8*
Overall Width in.	17 1/4*	17 1/4*	17 3/4*	17 3/4*	23 3/4*	23 3/4*	23 3/4*	23 3/4*
Overall Depth in.	19 1/4*	19 1/4*	19 1/2*	19 1/2*	24 1/4*	24 1/2*	25 1/2*	24 1/2*
Interior cu. ft.	1	1	2.5	2.5	5.2	4.6	4.6	4.5
Number of Shelves	2, adj	2 adj	4 adj	4 adj	3 adj	2 adj, 1 fixed	2 adj, 1 fixed	2 adj, 1 fixed
Door Type	solid	glass	solid	glass	glass	solid	glass	Solid SS
Defrost Type	cycle	cycle	cycle	cycle	cycle	cycle	cycle	cycle
Compressor Size	1/8 HP	1/8 HP	1/8 HP	1/8 HP	1/6HP	1/8 HP	1/8 HP	1/8 HP
Plug (Receptacle)	5-15P	5-15P	5-15P	5-15P	5-15P	5-15P	5-15P	5-15P
Total Amp Draw	0.9	0.9	0.9	0.9	1.3	1.74	1.74	1.74
Voltage	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz

Left hinged models available!

**ADA model available

*Exterior dimensions do not include 4" clearance required for installation.
 Specifications subject to change without notice.

Receptacle Required



115V/1PH/ 60Hz
 NEMA-5-15P



LABRepCo

One Company • Many Solutions

Revision Date: 4/25/2016