



Silver Series Manual Defrost -45°C Chest Freezers

The LabRepCo FUTURA Silver Series manual defrost -45 °C chest freezers are designed for laboratory & pharmaceutical applications that require more rigorous temperature control and construction features. These durable, easy to use units feature an intelligent microprocessor controller with digital display and adjustable temperature control.

Standard Product Features:

- Digital temperature display
- Intelligent microprocessor controller with digital display and adjustable temperature control
- Visual and audible hi or low temperature, sensor failure, door open alarms
- Keyed door locks
- Magnetically sealed door with rapid cool-down mode after door opening for quick temperature recovery
- Cabinet constructed of white powder-coated 304 steel with CFC-free, high density insulation (2.4" thick)
- Easy swivel casters
- Doors: 1 Lid
- Defrost: Manual
- Built and shipped from ISO 9001:2008 manufacturing and distribution facilities
- Each cabinet undergoes three quality control checks prior to shipping
- One wire basket (L2X-7-C45 and L2X-11-C45 only)
- Certification: UL / cUL
- Warranty: Two-year parts & compressor



L2X-3-C45



L2X-7-C45



L2X-11-C45

Model	CAPACITY (cu. ft.)	Exterior Dimensions (HxWxD)	Interior Dimensions (HxWxD)	Temperature Range (Setpoint)	Shipping Weight
L2X-3-C45	3.5	33.5" x 27.25" x 22"	24.25" x 20.5" x 13.5"	-10 °C to -45°C °C	102 lbs
L2X-7-C45	7	35" x 37.13" x 30.5"	24" x 29.6" x 19.1"	Preset -40 °C	113.5 lbs
L2X-11-C45	11	33.5" x 50.25" x 30.5"	24" x 42.38" x 19"		176 lbs