

Futura
-20°C
-30°C
-40°C

Silver Series Undercounter Freezers

Our **FUTURA Silver Series** undercounter freezers are designed for applications in the research market requiring safe and reliable sample storage in a space-saving design.



LHP-5-UFMB

Standard Product Features:

- Auto or Manual defrost
- Microprocessor temperature controller*
- Digital temperature display*
- Audible and visual high and low temperature alarms*
- Remote alarm contacts*
- Adjustable operating temperature range: -15°C to -25°C
- Keyed door lock
- Magnetic door gasket for positive seal
- 2 probes; 1 in air and 1 in glycol bottle*
- Back wall evaporator cover*
- Forced draft circulation*
- Leveling legs
- Probe access port (3/8")
- Powder coat exterior
- High density urethane foam cabinet and door insulation
- HFC-free refrigerant (R600a), foam insulation and packaging (LH-2-FM features R134a refrigerant)
- Warranty: 2 years parts and labor (LHP) models
1 year parts and labor (LH) models

*featured only on LHP models

Optional Accessories:

- NIST Traceable Thermometers
- NIST Traceable Data Loggers

LABRepCo

One Company • Many Solutions

Auto Defrost Models



Manual Defrost Models



Hydrocarbon Refrigerant used in all Silver Series undercounter refrigerators

Model	LH-2-FM	LHP-2-FM	LHP-2-FA	LHP-5-UFMB**	LHP-5-UFMBSS	LHP-5-UFAB	LHP-5-UFM30	LHP-5-UFM40
Weight lbs	46	40	40	100	100	100	93	93
Overall Height in.	19 1/2*	20 3/8*	20 3/8*	33 3/8*	33 3/8*	33 1/8*	32 1/8*	32 1/8*
Overall Width in.	18 1/2*	23 3/4*	23 3/4*	23 3/4*	23 3/4*	23 3/4*	23 3/4*	23 3/4*
Overall Depth in.	19 1/2*	24*	24*	24 1/2*	24 1/2*	24*	24*	24*
Interior cu. ft.	1.5	1.7	1.7	4.2	4.5	4.2	4	4
Avg. Temperature	-20°C	-20°C	-20°C	-20°C	-20°C	-20°C	-30°C	-40°C
# of Shelves	1, adj	2, adj	2, adj	2, fixed	2, fixed	3, adj	3, adj	3, adj
Door Type	solid	solid	solid	solid	solid	solid	solid	solid
Defrost Type	Manual	Manual	Auto	Manual	Manual	Auto	Manual	Manual
Compressor Size	1/10 HP	1/8 HP	1/8 HP	1/8 HP	1/8 HP	1/8 HP	1/4HP	1/4 HP
Plug (Receptacle)	5-15P	5-15P	5-15P	5-15P	5-15P	5-15P	5-15P	5-15P
Total Amp Draw	1.1	1	1	1.5	1.5	1.8	2.2	2.2
Voltage	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz

Please Note: All specifications are subject to change without notice.

**ADA model available

*Exterior dimensions do not include 4" clearance required for installation.

Specifications subject to change without notice.

Receptacle Required



115V/1PH/ 60Hz
NEMA-5-15P



LABRepCo

One Company • Many Solutions

Revision Date: 4/25/2016