



Stainless Steel Upright Laboratory Refrigerators

The **Futura Silver Series PRIME** stainless steel refrigerators are designed for a variety of applications in the research market requiring more rigorous temperature control and construction features. These units feature all stainless steel exteriors and interior floors, microprocessor temperature controller with audible and visual alarms, remote alarm contact, digital temperature display, key door locks and swivel casters (field installed).



LHP-23-HGSS

Standard Product Features:

- 18 & 30 Cu. Ft. Capacities
- Cycle defrost
- Adjustable mechanical thermostat
- 1°C to 10°C operating temperature
- Magnetic door gasket for positive seal
- Stainless steel exterior
- Swinging door(s) with heavy duty key lock
- Swivel casters
- Audible and visual Hi/Lo temperature alarm
- Digital temperature display
- Hydrocarbon, natural refrigerants (R290)
- Forced draft air circulation
- LED Interior lights
- 115V, 60Hz operation
- Probe access port (3/8")
- ETL Listed
- Warranty: 2 years parts & labor, 3 years compressor parts only

Optional Features:

- Chart Recorder
- Duplex Outlet
- Additional shelves
- 2" Access Ports



Specifications and Ordering Information:

Models	LHP-23-SDSS	LHP-23-HGSS	LHP-49-SDSS	LHP-49-HGSS	LHP-72-SDSS	LHP-72-HGSS
Weight (lbs.)	295	356	450	557	666	786
Overall Height (in.)	83 7/8	83 1/4	84	83 1/4	83 1/4	83 1/4
Overall Width (in.)	26 3/4	26 3/4	54	54	81	81
Overall Depth (in.)	34 7/8	35 3/4	35	35 3/4	35 3/4	35 3/4
Interior Cu. Ft.	23	23	49	49	72	72
Number of Shelves	3 adj	4 adj	6 adj	8 adj	9 adj	12 adj
Defrost Type	Cycle	Cycle	Cycle	Cycle	Cycle	Cycle
Refrigerant	R290a	R290a	R290a	R290a	R290a	R290a
Compressor Size HP	1/5	1/5	1/4	3/4	1/2	1/2
Plug (receptacle)	5-15R	5-15R	5-15R	5-15R	5-15R	5-15R
Total Amp Draw	3	3	4.5	4.5	7.5	7.5
Voltage	115V, 1PH, 60Hz	115V, 1PH, 60Hz	115V, 1PH, 60Hz	115V, 1PH, 60Hz	115V, 1PH, 60Hz	115V, 1PH, 60Hz

Please Note: Interior dimensions are estimates.
All specifications subject to change without notice.

LHP-72-SDSS



Receptacle Required



115V/1PH/ 60Hz
NEMA-5-15P



LABRepCo