



**Labsit Chair Plus 9123+,
Desk-to-Bench Height,
Freewheeling Casters,
Black Vinyl, Foot Ring,
With Arms**

Article #: LRC-LBS23PA*

**In-Stock Quick Ship Item*

Features

Labsit is a progressive chair design for all types of laboratories. The intelligent inherent flex of the seat shell offers comfort and ergonomic support with minimum adjustment. Lightweight and easy to maneuver, Labsit is a winning solution for lab stations, ESD areas and classrooms with an unbeatable value and performance.

Specifications

Seat Adjustment:	Desk-to-Bench Height Range 19-2/3" - 27-1/2"
Seat Depth:	18-1/8"
Shell Colors:	In-Stock (Blue); Factory Orders (Orange, Green, White, Black)
Seat Material:	Black Vinyl
Castors:	Freewheeling
Base:	Black Plastic Base
Cleanroom Rating:	ISO Level 5
Bio-Safety BSL:	Level 2
BIFMA Rating:	Complies with Level e3 Platinum
SEFA Rating:	Complies with Scientific Equipment & Furniture Assn. - SEFA12
Weight Capacity:	300 lbs.
Product Weight:	25 lbs.
Warranty:	10-Year Warranty
Original Design:	German Engineering and Production
ESD Capable:	Option Available by Request