



**Neon Series Chair Plus
9573+, Desk-to-Bench
Height, Freewheeling
Casters, Foot Ring,
Black Vinyl, Armless**

Article #: LRC-NEO73P*

**In-Stock Quick Ship Item*

Features

- The award-winning Neon Chair is known as the best lab chair globally for its design, superior quality, ergonomic features and overall comfort.
- Innovative Seat Cushion System allows for easy removal and replacements.
- All-inclusive ergonomic functionality with premium mechanism, back height lumbar adjustment, back lock, seat height adjustment, swivel-tilt, forward tilt and tilt tension.

Specifications

Seat Adjustment:	Desk-to-Bench Height Range 21-2/3" - 29-1/8"
Seat Depth Range:	16-5/8" - 18-7/8"
Accent Trim Color:	In-Stock (Blue); Factory Orders (Orange, Green and Grey)
Frame Color:	Grey Backs and Blue Trim
Seat Material:	Black Vinyl and Enhanced Comfort Seat
Castors:	Freewheeling
Base:	Black Metal Base
Cleanroom Rating:	ISO Level 7
Bio-Safety BSL:	Level 2
BIFMA Rating:	Complies with Level e3 Platinum
Lab Grade Rating:	Complies with Scientific Equipment & Furniture Assn. - SEFA12
Weight Capacity:	300 lbs.
Product Weight:	40 lbs.
Warranty:	10-Year Warranty
Original Design:	German Engineering and Production
ESD Capable:	Available by Request