




Upright Freezer Racks













For 3.75" High Standard Boxes



Polycarbonate Plastic Box
(Check Page 110 for Boxes)


- ◆ Racks are for 3.75" high standard polycarbonate boxes.
- ◆ Made of corrosion resistant stainless steel.



Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Shipping Weight lb(Kg)
UF-22-375	Rack Only	11 ³ / ₈ x 8 ¹ / ₁₆ x 5 ¹ / ₂ inch 289 x 205 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	2 x 2 =4 	2.9(1.33)
UF-23-375	Rack Only	11 ³ / ₈ x 12 ¹ / ₁₆ x 5 ¹ / ₂ inch 289 x 307 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	2 x 3 =6 	3.9(1.78)
UF-24-375	Rack Only	11 ³ / ₈ x 16 ¹ / ₁₆ x 5 ¹ / ₂ inch 289 x 408 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	2 x 4 =8 	4.9(2.24)
UF-32-375	Rack Only	16 ¹³ / ₁₆ x 8 ¹ / ₁₆ x 5 ¹ / ₂ inch 427 x 205 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	3 x 2 =6 	4.1(1.87)
UF-33-375	Rack Only	16 ¹³ / ₁₆ x 12 ¹ / ₁₆ x 5 ¹ / ₂ inch 427 x 307 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	3 x 3 =9 	5.5(2.51)
UF-34-375	Rack Only	16 ¹³ / ₁₆ x 16 ¹ / ₁₆ x 5 ¹ / ₂ inch 427 x 408 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	3 x 4 =12 	6.9(3.15)
UF-42-375	Rack Only	22 ¹ / ₈ x 8 ¹ / ₁₆ x 5 ¹ / ₂ inch 562 x 205 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	4 x 2 =8 	5.3(2.42)
UF-43-375	Rack Only	22 ¹ / ₈ x 12 ¹ / ₁₆ x 5 ¹ / ₂ inch 562 x 307 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	4 x 3 =12 	7.1(3.24)
UF-44-375	Rack Only	22 ¹ / ₈ x 16 ¹ / ₁₆ x 5 ¹ / ₂ inch 562 x 408 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 136 x 100 x 137 mm	4 x 4 =16 	9.0(4.06)
UF-52-375	Rack Only	27 ³ / ₁₆ x 8 ¹ / ₁₆ x 5 ¹ / ₂ inch 691 x 205 x 140 mm	5 ¹ / ₄ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 133.5 x 100 x 137 mm	5 x 2 =10 	6.5(2.94)
UF-53-375	Rack Only	27 ³ / ₁₆ x 12 ¹ / ₁₆ x 5 ¹ / ₂ inch 691 x 307 x 140 mm	5 ¹ / ₄ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 133.5 x 100 x 137 mm	5 x 3 =15 	8.7(3.94)
UF-54-375	Rack Only	27 ³ / ₁₆ x 16 ¹ / ₁₆ x 5 ¹ / ₂ inch 691 x 408 x 140 mm	5 ¹ / ₄ x 3 ¹⁵ / ₁₆ x 5 ³ / ₈ inch 133.5 x 100 x 137 mm	5 x 4 =20 	10.9(4.94)


Upright Freezer Drawer Racks













For 3.75" High Standard Boxes



Polycarbonate Plastic Box
(Check Page 110 for Boxes)

- ◆ Racks are for 3.75" high standard polycarbonate boxes.
- ◆ These drawer racks facilitate both access to and storage of samples without having to remove the whole rack from freezer.
- ◆ Made of corrosion resistant stainless steel.



Part Number	Description	Dimension L x H x W	Compartment Size l x h x w	Capacity (Box)	Shipping Weight lb(Kg)
UFD-22-375	Rack Only	11 ¹ / ₄ x 8 ³ / ₈ x 5 ¹ / ₂ inch 286 x 213 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	2 x 2 =4	 4.3(1.97)
UFD-23-375	Rack Only	11 ¹ / ₄ x 12 ¹ / ₂ x 5 ¹ / ₂ inch 286 x 318 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	2 x 3 =6	 7.1(2.81)
UFD-24-375	Rack Only	11 ¹ / ₄ x 16 ¹¹ / ₁₆ x 5 ¹ / ₂ inch 286 x 424 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	2 x 4 =8	 8.1(3.66)
UFD-32-375	Rack Only	16 ⁵ / ₈ x 8 ³ / ₈ x 5 ¹ / ₂ inch 423 x 213 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	3 x 2 =6	 6.0(2.65)
UFD-33-375	Rack Only	16 ⁵ / ₈ x 12 ¹ / ₂ x 5 ¹ / ₂ inch 423 x 318 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	3 x 3 =9	 8.3(3.78)
UFD-34-375	Rack Only	16 ⁵ / ₈ x 16 ¹¹ / ₁₆ x 5 ¹ / ₂ inch 423 x 424 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	3 x 4 =12	 10.8(4.92)
UFD-42-375	Rack Only	22 ³ / ₁₆ x 8 ³ / ₈ x 5 ¹ / ₂ inch 563 x 213 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	4 x 2 =8	 7.4(3.34)
UFD-43-375	Rack Only	22 ³ / ₁₆ x 12 ¹ / ₂ x 5 ¹ / ₂ inch 563 x 318 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	4 x 3 =12	 10.5(4.76)
UFD-44-375	Rack Only	22 ³ / ₁₆ x 16 ¹¹ / ₁₆ x 5 ¹ / ₂ inch 563 x 424 x 140 mm	5 ³ / ₈ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 137 x 100 x 133 mm	4 x 4 =16	 13.6(6.19)
UFD-52-375	Rack Only	26 ⁷ / ₈ x 8 ³ / ₈ x 5 ¹ / ₂ inch 683 x 213 x 140 mm	5 ⁵ / ₁₆ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 134.5 x 100 x 133 mm	5 x 2 =10	 8.8(3.97)
UFD-53-375	Rack Only	26 ⁷ / ₈ x 12 ¹ / ₂ x 5 ¹ / ₂ inch 683 x 318 x 140 mm	5 ⁵ / ₁₆ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 134.5 x 100 x 133 mm	5 x 3 =15	 11.9(5.40)
UFD-54-375	Rack Only	26 ⁷ / ₈ x 16 ¹¹ / ₁₆ x 5 ¹ / ₂ inch 683 x 424 x 140 mm	5 ⁵ / ₁₆ x 3 ¹⁵ / ₁₆ x 5 ¹ / ₄ inch 134.5 x 100 x 133 mm	5 x 4 =20	 16.2(7.34)