

# myFUGE™ 5



- + For quick spin downs of 5ml centrifuge tubes
- + Includes adapters for 12 x 75mm tubes
- + Starts & stops with closing/opening of the lid
- + Conserves valuable bench space

CE



Combination rotor accepts  
1.5/2.0ml & 5ml centrifuge tubes



New 5ml Centrifuge Tubes  
(Item: C1005-T5)

The myFuge™ 5 microcentrifuge includes a four position rotor for compatibility with a variety of 5ml test tubes, including 5ml snap cap centrifuge tubes. In addition, this combination rotor has four additional slots for centrifuging standard microtubes (1.5/2.0ml). Simply close the lid and the centrifuge quickly accelerates to the fixed speed of 5,500rpm. Open the lid, and the rotor decelerates to a stop for removing samples.

The 5ml rotor positions accepts 5ml centrifuge tubes (snap cap) directly. While 5ml blood and culture tubes (12x75mm) and cryovials are also accepted with the included adapters. The microtube positions are designed to accept 0.5 and 0.2ml tubes with the use of optional adapters.

## Technical Data:

|             |                                        |
|-------------|----------------------------------------|
| Speed:      | 5,500rpm / 2,000 xg                    |
| Capacity:   | 4 x 5ml tubes & 4 x 1.5/2.0ml          |
| Dimensions: | 5.5 x 7.9 x 5.5 in.<br>14 x 20 x 14 cm |
| Weight:     | 2.2 lbs. / 1 kg                        |
| Electrical: | 115V 50Hz or 230V, 60Hz                |
| Warranty:   | 1 Year                                 |



## Ordering Information:

|               |                                                                                     |
|---------------|-------------------------------------------------------------------------------------|
| <b>C1005*</b> | <b>myFuge 5 Mini Centrifuge with combination rotor, 4 x 5ml &amp; 4 x 1.5/2.0ml</b> |
| C1008-A5-6    | 0.5ml Adapter Pack, 6/pk                                                            |
| C1008-A2-6    | 0.2ml Adapter Pack, 6/pk                                                            |
| C1005-AC2     | Adapter Pack, for Cryovials (0.5 to 2.0ml) and 12x32mm HPLC vials, 4/pk             |
| C1005-T5      | 5ml centrifuge tube, clear, 200/pk (2 bags of 100)                                  |
| C1005-T5-ST   | 5ml <i>sterile</i> centrifuge tube, clear, 200/pk. (10 bags of 20)                  |

\*115V with US plug. To order in 230V, add (-E) to the item number.

**Benchmark**  
Scientific

**Benchmark Scientific Inc.**  
PO Box 709 Edison, NJ 08818  
Phone: 908-769-5555 Fax: 908-222-1864  
Web: [www.BenchmarkScientific.com](http://www.BenchmarkScientific.com)  
Email: [Info@BenchmarkScientific.com](mailto:Info@BenchmarkScientific.com)