

The CliniCool® Silver Series PRIME stainless steel pharmacy refrigerators are designed for applications in the clinical, pharmacy, and medical market requiring more rigorous temperature control and construction features. These units features solid stainless steel doors with pull handles, all stainless steel exterior, microprocessor temperature controller with audible and visual alarms, remote alarm contacts, digital temperature display, keyed door locks, swivel casters (field installed), pharmacy toolkit and CDC approved power cords.



LHP-23-SDSS-PH

Standard Product Features:

- Cycle defrost
- Microprocessor temperature controller
- Adjustable operating temperature: 36° F to 46° F
- Digital Data Logger (DDL) meets current CDC requirements for vaccine storage and monitoring. F/C switchable, has a 3-year certification of calibration, and (1) buffered probe in product simulated solution. Min/Max memory along with Alarm event handling. USB port for downloading stored data
- Audible and visual temperature alarms
- Remote alarm contacts
- Magnetic door gasket for positive seal
- Swivel Casters
- White cabinet and door
- LED Interior Lights
- High density urethane foam cabinet and door insulation
- Hydrocarbon, natural refrigerant (R290)
- 115V, 60Hz operation
- Probe access port (3/8")
- ETL Listed
- Warranty: 2 years parts & labor, 3 years compressor parts only



Hydrocarbon refrigerant used in all CliniCool Stainless Steel models

Vaccine Compliant Features

All units tested and certified to maintain temperatures within refrigerated or frozen ranges required by the CDC

(1) Data Logger with 1 bottle probe, 3 years NIST Certified of Calibration

Vaccine Storage Power Cord Warning Label

Vaccine Toolkit

Specifications and Ordering Information:

Models	LHP-23-SDSS-PH	LHP-23-HGSS-PH	LHP-49-SDSS-PH	LHP-49-HGSS-PH	LHP-72-SDSS-PH	LHP-72-HGSS-PH
Weight (lbs.)	295	356	450	557	666	786
Overall Height (in.)	83 7/8	83 1/4	84	83 1/4	83 1/4	83 1/4
Overall Width (in.)	26 3/4	26 7/8	54	54	81	81
Overall Depth (in.)	34 7/8	35 3/4	35	35 3/4	35 3/4	35 3/4
Interior Cu. Ft.	23	23	49	49	72	72
Number of Shelves	3 adj	3 adj	6 adj	6 adj	9 adj	9 adj
Defrost Type	Cycle	Cycle	Cycle	Cycle	Cycle	Cycle
Refrigerant	R290a	R290a	R290a	R290a	R290a	R290a
Compressor Size HP	1/5	1/5	1/4	1/4	1/2	1/2
Plug (receptacle)	5-15R	5-15R	5-15R	5-15R	5-15R	5-15R
Total Amp Draw	3	3	4.5	4.5	7.5	7.5
Voltage	115V, 1PH, 60Hz	115V, 1PH, 60Hz	115V, 1PH, 60Hz	115V, 1PH, 60Hz	115V, 1PH, 60Hz	115V, 1PH, 60Hz

Please Note: Interior dimensions are estimates.
All specifications subject to change without notice.

LHP-72-HGSS-PH



Receptacle Required



115V/1PH/ 60Hz
NEMA-5-15P



LABRepCo
An Employee-Owned Company