



Exterior General Purpose Refrigerator LABH-30-SD

Front View

Side View

Top View



115VAC/60Hz/1PH
NEMA-5-15R

*Leveling legs are fully recessed

Model #	Cu. Ft	Defrost	Door	Int Door	Shelves	W"	D"	H"	H.P.	Amps	Weight
LABH-30-SD	30	Auto	1 Solid	0	5 Adj	$35\frac{3}{4}$	$34\frac{1}{4}$	$72\frac{1}{2}$ *	$\frac{1}{3}$	6.8	312 lbs.