



Pass-Thru Laboratory Refrigerators

The **LabRepCo FUTURA PLUS+ Series pass-thru refrigerators** are designed for many different applications in the research and clinical market requiring durable refrigeration plus convenient access of samples from side-by-side rooms. These refrigerators are available in 25, 50 and 80 cu. ft. sizes and come standard with glass doors in the front and back of the unit, and a stainless steel interior and exterior.



LABL-25-HGPT

Standard Product Features:

- Cycle, auto defrost
- +1 °C to +8 °C average operating temperature
- Temperature controller with digital display and Hi/Low alarm
- CFC-free refrigerant (R134a), polyurethane foamed in place insulation
- Heavy duty epoxy coated, wire shelves
- Interior fluorescent lights
- Fin and tube, forced air evaporator
- Forced air circulation
- Fin and tube compressor
- Adjustable thermostat
- All stainless steel interior and exterior with glass door(s) in front and back
- 4 heavy duty 3" swivel casters (locking)
- Key lock doors
- 115V operation
- Warranties: 3 years parts and labor, additional 2 years compressor parts only



**Microprocessor Temperature
Controllers**





LABL-50-HGPT



LABL-80-HGPT

Models	LABL-25-HGPT	LABL-50-HGPT	LABL-80-HGPT
Crated Weight lbs.	500	650	850
Overall Height in.	81 1/2	81 1/2	81 1/2
Overall Width in.	26	52	78
Overall Depth in.	37 1/2	37 1/2	37 1/2
Interior Height in.	59	59	59
Interior Width in.	23	48	75
Interior Depth in.	30	30	30
Interior cu. ft.	25	50	80
Number of Doors	2	4	6
Number of Shelves	3	6	9
Number of Casters	4	4	4
Compressor Size hp	1/3	3/4	1
Refrigerant	R134a	R134a	R134a
Running Amps/Breaker	6.0/15	12.0/15	14.0/15
Voltage Description	115V, 1PH, 60 Hz	115, 1PH, 60 Hz	115, 1PH, 60 Hz
Plug Configuration	5-15P NEMA	5-15P NEMA	5-15P NEMA

Please Note: All dimensions are nominal. All specifications subject to change without notice.

Receptacle Required



115V/1PH/ 60Hz
NEMA-5-15R

LABRepCo
An Employee-Owned Company