



MetroMax iQ™ Storage System

MetroMax Q™ Shelving Starter and Add-on Units

- Save time and money by eliminating adjacent posts and replacing with "S" hooks.
- Starter Units consist of four posts and indicated number of MetroMax Q shelves.
- Add-On Units consist of two posts, indicated number of shelves, and two "S" Hooks per shelf, for attaching shelf to adjacent unit. To attach shelves at right angles, Intermediate "S" Hook Kits must be ordered separately.



Add-on units with "S" hooks can only be attached to shelves supported by two or more posts.

"S" Hook

Add on shelving using two posts instead of four. Includes one "S" hook and one corner collar plug. Two kits required per level when configuring units end-to-end.



"S" Hook

MetroMax Q™ Intermediate "S" Hook Kit

Use when configuring MetroMax Q Starter and Add-On Units at right angles.



Intermediate "S" Hook Kit

Shelving Starter and Add-On Units — Easy Order Guide

5 Shelves
74" (1880mm) Posts

	Cat. No.	18" (457mm) Wide		21" (530mm) Wide		24" (610mm) Wide	
		Starter	Add-On	Starter	Add-On	Starter	Add-On
24" (610mm)		LAB3G713Q5	LAB3G713QA5	LAB3G714Q5	LAB3G714QA5	LAB3G715Q5	LAB3G715QA5
Long	Wt. lbs. (kg)	50.4 (23)	40.4 (18)	56.4 (26)	46.4 (21)	56.4 (26)	46.4 (21)
30" (760mm)		LAB3G723Q5	LAB3G723QA5	LAB3G724Q5	LAB3G724QA5	LAB3G725Q5	LAB3G725QA5
Long	Wt. lbs. (kg)	57.9 (26)	47.9 (22)	63.1 (29)	53.1 (24)	66.4 (30)	56.4 (26)
36" (914mm)		LAB3G733Q5	LAB3G733QA5	LAB3G734Q5	LAB3G734QA5	LAB3G735Q5	LAB3G735QA5
Long	Wt. lbs. (kg)	65.4 (30)	55.4 (25)	70.9 (32)	60.9 (28)	76.4 (35)	65.4 (30)
42" (1066mm)		LAB3G743Q5	LAB3G743QA5	LAB3G744Q5	LAB3G744QA5	LAB3G745Q5	LAB3G745QA5
Long	Wt. lbs. (kg)	72.9 (33)	62.9 (29)	82.1 (37)	72.1 (33)	86.4 (39)	78.4 (36)
48" (1219mm)		LAB3G753Q5	LAB3G753QA5	LAB3G754Q5	LAB3G754QA5	LAB3G755Q5	LAB3G755QA5
Long	Wt. lbs. (kg)	80.4 (36)	70.4 (32)	89.4 (41)	79.4 (36)	96.4 (44)	86.4 (39)
60" (1524mm)		LAB3G763Q5	LAB3G763QA5	LAB3G764Q5	LAB3G764QA5	LAB3G765Q5	LAB3G765QA5
Long	Wt. lbs. (kg)	95.4 (43)	85.4 (39)	106.4 (48)	96.4 (44)	116.4 (53)	106.4 (48)
72" (1825mm)		LAB3G773Q5	LAB3G773QA5	LAB3G774Q5	LAB3G774QA5	LAB3G775Q5	LAB3G775QA5
Long	Wt. lbs. (kg)	110.4 (50)	100.4 (46)	121.4 (55)	111.4 (51)	136.4 (62)	126.4 (57)