

Worktables — 17.20

Built for cleanliness and stability.

- Seamless 14-gauge Type 304 stainless steel work surface and support structure.
- Shipped knocked down and can be assembled in minutes without tools.

All Worktables are 34" (864mm) high with stationary posts and leveling feet.

Standard Worktables — 30" (760mm) Wide

Worktable with Bottom Shelf

Length (in.) (mm)	Post Material	Bottom Shelf Material	Bottom 3-Sided Frame Material	Bottom H-Frame Material	Approx. Pkd. Wt. (lbs.) (kg)	Cat. No.
48 1219	SS	SS	—	—	90 41	LABSF503TW
60 1524	SS	SS	—	—	130 59	LABSF603TW
72 1825	SS	SS	—	—	170 77	LABSF703TW
96 2439	SS	SS	—	—	226 103	LABSF903TW
60 1524	Chrome	Galvanized	—	—	146 66	LABCF603TW
72 1825	Chrome	Galvanized	—	—	180 82	LABCF703TW
96 2439	Chrome	Galvanized	—	—	226 103	LABCF903TW
48 1219	SS	—	SS	—	115 52	LABSU503TW
60 1524	SS	—	SS	—	115 52	LABSU603TW
72 1825	SS	—	SS	—	142 65	LABSU703TW
96 2439	SS	—	SS	—	201 91	LABSU903TW
60 1524	SS	—	—	SS	115 52	LABSH603TW
72 1825	SS	—	—	SS	142 65	LABSH703TW
96 2439	SS	—	—	SS	202 92	LABSH903TW

SS = Stainless Steel



Worktable with 3-Sided Frame

Standard Worktables — 36" (914mm) Wide

Length (in.) (mm)	Post Material	Bottom Shelf Material	Bottom 3-Sided Frame Material	Bottom H-Frame Material	Approx. Pkd. Wt. (lbs.) (kg)	Cat. No.
60 1524	SS	SS	—	—	125 57	LABSF663TW
72 1825	SS	SS	—	—	130 59	LABSF763TW
96 2439	SS	SS	—	—	185 84	LABSF963TW
60 1524	Chrome	Galvanized	—	—	125 57	LABCF663FC
72 1825	Chrome	Galvanized	—	—	130 59	LABCF763FC
96 2439	Chrome	Galvanized	—	—	185 84	LABCF963FC
60 1524	SS	—	SS	—	104 47	LABSU663TW
72 1825	SS	—	SS	—	106 48	LABSU763TW
96 2439	SS	—	SS	—	161 73	LABSU963TW
60 1524	SS	—	—	SS	103 47	LABSH663TW
72 1825	SS	—	—	SS	107 49	LABSH763TW
96 2439	SS	—	—	SS	162 74	LABSH963TW

SS = Stainless Steel



Worktable with H-Frame

Standard Worktables — 44" (1118mm) Wide

Length (in.) (mm)	Post Material	Bottom Shelf Material	Bottom 3-Sided Frame Material	Bottom H-Frame Material	Approx. Pkd. Wt. (lbs.) (kg)	Cat. No.
60 1524	SS	SS	—	—	136 62	LABSF644TW
72 1825	SS	SS	—	—	160 73	LABSF744TW
96 2439	SS	SS	—	—	202 92	LABSF944TW
60 1524	Chrome	Galvanized	—	—	136 62	LABCF644TW
72 1825	Chrome	Galvanized	—	—	160 73	LABCF744TW
96 2439	Chrome	Galvanized	—	—	202 92	LABCF944TW
60 1524	SS	—	SS	—	114 52	LABSU644TW
72 1825	SS	—	SS	—	135 61	LABSU744TW
96 2439	SS	—	SS	—	177 81	LABSU944TW
60 1524	SS	—	—	SS	113 51	LABSH644TW
72 1825	SS	—	—	SS	136 62	LABSH744TW
96 2439	SS	—	—	SS	178 81	LABSH944TW

SS = Stainless Steel

Mobile Worktables — 17.20

All Mobile Worktables are 34" (864mm) high. Mobile units have special posts and four 5" (127mm) poly casters — 2 swivel/2 brake.

Mobile Worktables — 30" (760mm) Wide

Length (in.) (mm)	Post Material	Bottom Shelf Material	Bottom 3-Sided Frame Material	Bottom H-Frame Material	Approx. Pkd. Wt. (lbs.) (kg)	Cat. No.
48 1219	SS	SS	—	—	105 48	LABSF503TWM
60 1524	SS	SS	—	—	160 73	LABSF603TWM
72 1825	SS	SS	—	—	182 83	LABSF703TWM
96 2439	SS	SS	—	—	226 103	LABSF903TWM
60 1524	Chrome	Galvanized	—	—	146 66	LABCF603TWM
72 1825	Chrome	Galvanized	—	—	180 82	LABCF703TWM
96 2439	Chrome	Galvanized	—	—	226 103	LABCF903TWM
48 1219	SS	—	SS	—	130 59	LABSU503TWM
60 1524	SS	—	SS	—	139 63	LABSU603TWM
72 1825	SS	—	SS	—	155 71	LABSU703TWM
96 2439	SS	—	SS	—	201 91	LABSU903TWM
60 1524	SS	—	—	SS	124 56	LABSH603TWM
72 1825	SS	—	—	SS	156 71	LABSH703TWM
96 2439	SS	—	—	SS	202 92	LABSH903TWM

SS = Stainless Steel



Mobile Worktable

Worktables with Overhead — 17.26

All Worktables with Overhead are 30" (760mm) wide and have rear posts that are approximately 88" (2235mm) high. Cantilever shelves for overhead models are sold separately (see page 194).

Length (in.) (mm)	Post Material	Bottom Shelf Material	Bottom 3-Sided Frame Material	Approx. Pkd. Wt. (lbs.) (kg)	Cat. No.
60 1524	SS	SS	—	158 71	LABSF603CTW
72 1825	SS	SS	—	251 114	LABSF703CTW
96 2440	SS	SS	—	288 131	LABSF903CTW
60 1524	Chrome	Galvanized	—	214 97	LABCF603CTW
72 1825	Chrome	Galvanized	—	251 114	LABCF703CTW
96 2440	Chrome	Galvanized	—	288 131	LABCF903CTW
60 1524	SS	—	SS	196 89	LABSU603CTW
72 1825	SS	—	SS	230 105	LABSU703CTW
96 2440	SS	—	SS	267 121	LABSU903CTW

All overhead models are fitted with stationary posts and leveling feet. For overhead models in other sizes, contact your InterMetro representative.

SS = Stainless Steel



Shown with optional Cantilever Shelf and Utility Rack