



## MultiWash+ Microplate Washer

A 96- and 384-well automated washer

### Benefits

- Walk-up usability
- Start simply with wash and rinse bottles (included)
- Flexible system washes 96- or 384-well plates
- Automated rinse reduces clogging
- Quiet and efficient

MultiWash™+ is an automated microplate washer that provides maximum flexibility with a minimum footprint. The MultiWash+ Microplate Washer can be configured for both 96- and 384-well plates and comes with four wash/rinse bottles and one waste bottle for out-of-the-box use.

### Wash Protocols & Modes

The on-board touch panel comes with the option of 20 different wash protocols, including up to 8 cycles within each protocol for walk-up usability and flexibility. Washing variations include adjustable speed and volume, adjustable aspiration speed and time, adjustable soak times and three modes of shaking are available to mix solutions.

### Efficiency

Efficient washing is achieved with cross-wise aspiration-reducing residual volume within each well. The washer has a small footprint and is vacuum- and pressure free, with on-board pumps that provide a quiet wash experience.

## Ordering information

- MultiWash+ for 96-well plates: 4 liquid lines, 8 way manifold and tubing, 4 wash/rinse bottles (2L), 1 waste bottle (2.5L)  
Part Number: MultiWash+
- Optional 16-way manifold  
Part Number: 5032338
- Optional 12-way manifold  
Part Number: 5032336
- Washer Maintenance Kit: 2 X T 2.5A spare fuses, 2 cleaning needles, tubing sets for all bottles, cleaning tool for needles, tubing set for manifold, and spare tubing set for manifold  
Part Number: 5032339
- External Tubing Kit (4-way)  
Part Number: 5032340
- Wash bottle (2L)  
Part Number: 5032341
- Rinse bottle (2L)  
Part Number: 5032342

## Technical Specifications

Dispense volume	50 – 200µL in 50µL increments
Precision	<5% at 300µL across the plate
Residual volume	<1µL per well
User interface	2 X 16 characters backlit display keyboard with 5 function keys
Shaking	3 speeds
Microplates	96-well plates flat and round bottom (384-well plates with 16-manifold) Fits plates with standard ANSI SBS footprint
Power requirement	90 – 250B, 60V auto-sensing, 50/60 Hz
Instrument dimensions (WxHxD)	21 x 21 x 44cm, 8.3" x 8.3" x 17.3"
Weight	8Kg (17.8lbs)

### Contact Us

Phone: +1-800-635-5577  
Web: [www.moleculardevices.com](http://www.moleculardevices.com)  
Email: [info@moldev.com](mailto:info@moldev.com)

Check our website for a current listing of worldwide distributors.

### Regional Offices

USA and Canada +1-800-635-5577  
Brazil +55-11-3616-6607  
China (Beijing) +86-10-6410-8669  
China (Shanghai) +86-21-3372-1088  
Germany 00800-665-32860

Japan (Osaka) +81-6-7174-8831  
Japan (Tokyo) +81-3-6362-5260  
South Korea +82-2-3471-9531  
United Kingdom +44-118-944-8000

