

CliniCool® Silver Series undercounter pharmacy/vaccine freezers are ideal for clinical labs, medical offices and pharmacies with a limited space and budget. These units offer superior temperature performance for safe storage of vaccines while meeting federal and state guidelines.


Standard Product Features:

- Auto or Manual defrost
- Microprocessor temperature controller with battery back-up
- Adjustable operating temperature range: +36°F to +46°F
- Freestanding thermometer with 2 year certificate of calibration
- Digital temperature display
- Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable
- Audible and visual high and low temperature alarms, with battery back up
- Remote alarm contacts
- 2 probes; 1 in air and 1 in glycol bottle
- Adjustable or fixed shelves with guard rail on back to prevent product from being pushed against the back wall
- Bottom footed shelf to allow air circulation
- Back wall evaporator cover
- Forced draft circulation
- Keyed door lock
- LED light
- Magnetic door gasket for positive seal
- Leveling legs
- Probe access port (3/8")
- Powder coat exterior
- High density urethane foam cabinet and door insulation
- HFC-free refrigerant (R600a), foam insulation and packaging
- Pharmacy refrigerator/freezer toolkit and temperature logs
- Pyxis®, Omnicell® and AcuDose RX® compatible
- Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage.
- Warranty: 2 years parts and labor



Hydrocarbon refrigerant used in all LHP medical freezer models

Vaccine Compliant Features




Specifications and Ordering Information:

Model	LHP-1-FA-PH	LHP-2-FM-PH*	LHP-5-UFMB-PH*	LHP-5-UFMBPH-ADA*	LHP-5-UFMBSS-PH*	LHP-5-UFAB-PH*
Weight lbs	80	80	100	108	100	100
Overall Height in.	20 3/8*	20 3/8*	33 3/8*	31 15/16*	33 3/8*	33 1/8*
Overall Width in.	23 3/4*	23 5/8*	23 3/4*	23 3/4*	23 3/4*	23 3/4*
Overall Depth in.	24*	24*	24 1/2*	24 1/2*	24 1/2*	24 1/2*
Interior Height in.	12 3/8	13 7/8	24 1/12	24	24	23 3/4
Interior Width in.	17 3/4	17 3/4	19	19	19	18 5/8
Interior Depth in.	7 1/2	16	17 1/4	17 1/4	17 1/4	16 1/4
Interior cu. ft.	1	1.7	4.2	4.2	4.2	4.2
Temperature Range	+5°F to -13°F	+5°F to -13°F	+5°F to -13°F	+5°F to -13°F	+5°F to -13°F	+5°F to -13°F
# of Shelves	2, adj	2, adj	2, fixed	2, fixed	2, fixed	3, adj
Door Type	solid	solid	solid	solid	solid	solid
Defrost Type	Control Auto	Manual	Manual	Manual	Manual	Auto
Compressor Size	1/8 HP	1/8 HP	1/8 HP	1/8 HP	1/8 HP	1/8 HP
Plug (Receptacle)	5-15P	5-15P	5-15P	5-15P	5-15P	5-15P
Total Amp Draw	1	1	1.5	1.5	1.5	1.5
Voltage	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz	115V, 60Hz

Please Note: All dimensions are nominal. All specifications subject to change without notice.

*Left hinge model available

Auto Defrost Models



Manual Defrost Models



Receptacle Required



15V/1PH/ 60Hz
IEMA-5-15P



LABRepCo
An Employee-Owned Company

101 Witmer Road, Suite 700 Horsham, PA 19044 | www.labrepco.com | (800) 521-0754 | sales@labrepco.com