



## -20°C Manual Defrost Laboratory Chest Freezers

The **FUTURA Silver Series** manual defrost freezers are designed for many different applications in the research and clinical market requiring more rigorous temperature control and construction features. They are available in several sizes that fit any storage requirements.

---

### **Standard Product Features:**

- Sizes from 15 to 20 cu. ft.
- Manual defrost
- Adjustable mechanical thermostat
- -20°C average operating temperature
- Magnetic door gasket for positive seal
- Solid door with key lock
- White cabinet and door
- Shelves: Two removable baskets
- High density urethane foam cabinet and door insulation
- CFC-Free refrigerant, foam insulation and packaging
- 115V, 60Hz operation
- Probe access port (3/8")
- CSA or UL/C-UL Listed
- Warranty: 2 years parts & labor, 5 years compressor parts only



**LABH-15-CFM**

### **Optional Features:**

- Casters\*
- Chart recorder
- Additional shelves\*\*

*\*Only available on 20 and 30 cu. ft. models*

*\*\*Available for LABH-30-FM model only*





LABH-15-CFM

Model	LABH-15-CFM	LABH-20-CFM
Net Weight lbs.	197	266
Overall Height in.	35 1/8	34
Overall Width in.	56	74 3/8
Overall Depth in.	29 1/2	29 1/2
Interior Height in.		24 1/2
Interior Width in.		67 1/2
Interior Depth in.		21 1/2
Interior cu. ft.	15	20
Temperature	-15°C to -25°C	-15°C to -25°C
Number of Shelves	N/A	N/A
Doors (Ext/Int)	1/0	1/0
Compressor Size HP	1/4	1/4
Receptacle	NEMA 5-15R	NEMA 5-15R
Total Amp Draw	5	5
Voltage	115V, 60Hz	115V, 60Hz
Listing	UL/C-UL	UL/C-UL

*Please Note: Interior dimensions are estimates.  
All specifications subject to change without notice.*



LABH-20-CFM

Receptacle Required



115V/1PH/ 60Hz  
NEMA-5-15P

