

# TW 34HC

## LN<sub>2</sub> HIGH-CAPACITY SMALL FREEZER

For long-term storage of large quantities at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Numbered index ring and tapered internal canister locator
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	34HC
LN <sub>2</sub> Capacity (L)	34
Static Hold Time (Days)	200
Evap. Rate (Liters/Day)	0.17
DIMENSIONS	
Necktube (in)	3.6
Overall Dia. (in)	18.8
Height (in)	26.3
WEIGHT	
Empty (Lbs)	35.38
Full Weight (Lbs)	95.9
CAPACITY	
Number of Canisters	6
Canister Dimensions (in)	2.75 x 11
1.2ml & 2.0 ml vials (5/cane)	570
1.2ml & 2.0 ml vials (6/cane)	684
.5 cc Straws (10/cane)	1920
.5 cc Straws- Bulk (2 Levels)	4800
.25 cc Straws (28/cane)	5376
.25 cc Straws- Bulk (2 Levels) in Goblets	10740



800.886.2796 | Fax 317. 297.7988  
[sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)  
[www.ICBiomedical.com](http://www.ICBiomedical.com)

© 2021, IC Biomedical, LLC