

Our FUTURA Silver Series dual temperature refrigerators/ freezers are designed for many different applications in the research and scientific markets. This unit features one solid exterior door with a stainless steel pull handle, one solid interior door, a mechanical thermostat and leveling legs.

SKU: LABH-30-DTM



Optional Features:

- Additional full shelf
- Chart recorder
- Casters
- Digital display

Specifications:

- Cu. Ft.: 30
 - Refrigerator: 25 Cu. Ft.
 - Freezer: 5 Cu. Ft.
- Door: single solid, white
- Temperature Range:
 - Refrigerator: 0 °C to +10 °C
 - Freezer: -15 °C to -20 °C
- Exterior Dimensions: 35.75"W x 35.5"D x 72.5"H
- Interior Dimensions:
 - Refrigerator: 32"W x 22.5"D x 50.25"H
 - Freezer: 29"W x 22.25"D x 12"H
- Alarm: N/A
- Casters: optional
- Shelves: 3 adjustable (refrigerator)
- Electrical Requirements: 115V, 60Hz, 1PH; 6 amps
- Supply Plug: 5-15P NEMA
- Shipping Weight: 383lbs

Standard Product Features:

- Manual defrost
- Adjustable manual thermostat
- One exterior & one interior door
- White powder coat paint cabinet and doors
- Magnetic door gasket for positive seal
- Heavy duty shelves
- Stainless steel pull handle
- High density urethane foam cabinet & door insulation
- CFC-free refrigerant, foam insulation & packaging
- Key door lock
- Solid doors, right hand swing
- Polyurethane foam insulation
- 115V operation
- Leveling Legs
- Two probe access ports (3/8"), one in each chamber
- Warranty: 2 years parts and labor, 5 year compressor parts only