

The LHP-11-HG-PH CliniCool Silver Series PRIME glass door pharmacy refrigerator is designed to securely store a wide range of vaccines & medications within strict temperature control between 36°F to 46°F (2°C to 8°C). This unit features a solid door refrigerator and comes standard with temperature alarms and is ENERGY Star Certified.

SKU: LHP-11-HG-PH



Optional Features:

- Additional full shelf
- NIST Traceable® Thermometers
- NIST Traceable® Data Loggers

Specifications:

- Cu. Ft.: 10.5
- Door: One glass swing door, right hinged
- Temperature Range:
+36°F to +46°F (+2°C to +8°C)
- Exterior Dimensions: 23.75"W x 26.5"D x 59"H
- Interior Dimensions: 19.87"W x 18.875"D x 53.5"H
- Alarm: Audible and visual high/low temperature alarms
- Casters: N/A
- Shelves: 6 adjustable
- Voltage: Domestic 115V, 60Hz
- Running Amps/Breaker: 2.3, 1/3HP
- Supply Plug: 5-15P NEMA
- Shipping Weight: 183 lbs.

Standard Product Features:

- ENERGY STAR Certified
- Cycle Defrost
- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- Adjustable operating temperature range:
+36°F to +46°F (+2°C to +8°C)
- Refrigerator/Freezer freestanding thermometer with 3 year NIST certificate of calibration, min/max memory, °F/°C switchable
- Shelves: 6 adjustable
- Bottom footed shelf with force air circulation
- Back wall evaporator cover
- Vaccine storage power cord warning label
- Pharmacy freezer tool kit & temperature logs
- Keyed door lock
- Magnetic door gasket for positive seal
- Leveling legs
- LED interior light
- Probe access port (3/8")
- Powder coated exterior and urethane foam insulation
- HFC-free refrigerant (R600a)
- Pyxis®, Omnicell® and AcuDose RX® compatible
- UL/C-UL Listed
- Warranty: 2 years parts and labor, additional 3 year compressor parts