

The LHP-11-SD FUTURA Silver Series compact solid door refrigerators are designed for a variety of research applications within the clinical market requiring strict temperature control. This unit features a solid door and comes standard with temperature alarms and is ENERGY STAR Certified.

SKU: LHP-11-SD



Optional Features:

- Additional full shelf
- NIST Traceable[®] Thermometers
- NIST Traceable[®] Data Loggers

Specifications:

- Cu. Ft.: 10.5
- Door: 1 solid swing door, right hinged
- Temperature Range:
+34 °F to +50 °F (+1 °C to +10 °C)
- Exterior Dimensions: 23.75"W x 26.5"D x 59"H
- Interior Dimensions: 19.87"W x 18.37"D x 53.5"H
- Alarm: Audible and visual high/low temperature alarms
- Casters: N/A
- Shelves: 6 adjustable
- Voltage: Domestic 115V, 60Hz, 1/3HP
- Running Amps/Breaker: 2.3
- Supply Plug: 5-15P NEMA
- Shipping Weight: 183 lbs

Standard Product Features:

- ENERGY STAR Certified
- Cycle Defrost
- Microprocessor temperature controller
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- Adjustable operating temperature range:
+34 °F to +50 °F (+1 °C to +10 °C)
- Shelves: 6 adjustable
- Keyed door lock
- Magnetic door gasket for positive seal
- Leveling legs
- LED interior lights, switch controlled
- Probe access port (3/8")
- Powder coat exterior
- High density urethane foam cabinet insulation
- HFC- free refrigerant (R600a)
- Warranty: 2 years parts and labor, additional 3 year compressor parts
- Forced draft air circulation
- UL/C-UL listed