



## Laboratory Chest Freezers

Our **FUTURA PLUS+ Series** manual defrost chest freezers are designed for applications in the standard research market requiring more rigorous temperature control and construction features.

---

### Standard Product Features:

- Manual defrost
- 5.0, 9.0 & 13.0 cu. ft. capacities
- 0°C to -25°C operating temperature
- Digital control with two temperature displays. One display shows the set temperature, the other shows the actual temperature in the freezer.
- Single-stage type refrigeration system with hermetic motor compressor
- Four heavy duty, swivel casters allow for convenient mobility
- White interior and exterior with hinged solid door
- 115V operation
- ETL electrical approved
- Warranties: 1 year parts and 90 days labor



LABL-5-CT25


### Optional Features:

- Chart recorder
- Alarm relay contacts



LABL-13-CT25


**LABRepCo**




**LABL-5-CT25**

Models	LABL-5-CT25	LABL-9-CT25	LABL-13-CT25
Overall Height in.	39	38 1/2	38 1/2
Overall Width in.	28 1/2	42	43
Overall Depth in.	22	24 1/2	31
Interior Height in.	27 1/2	28 1/2	28 1/2
Interior Width in.	23	36	37 1/2
Interior Depth in.	15	16 1/2	22
Interior cu. ft.	5.0	9.0	13.0
Number of Casters	4	4	4
Temperature	0° to -25°C	0° to -25°C	0° to -25°C
Voltage	115V, 1PH, 60Hz	115V, 1PH, 60Hz	115V, 1PH, 60Hz


Specifications subject to change without notice.



**LABL-9-CT25**



**LABL-13-CT25**



**Maximize Freezer Space**

Please visit [labrepco.com](http://labrepco.com) or speak with your local representative for inventory rack information.