

Heated and Cooled Incubators MIR and MLR Series

PHCbi brand heated and cooled incubators include a range of microprocessor based programmable temperature controlled cabinets. Choose from heated only, heated and refrigerated or a multi-parameter model with programmed lighting and controlled humidity.

■ Standard Feature □ Optional Feature

	Heated Incubator		Heated and Refrigerated Incubator			Heated, Refrigerated, Humidified, Lighted
	3.3 cu.ft.	5.4 cu.ft.	4.3 cu.ft.	8.4 cu.ft.	14.3 cu.ft.	10.4 cu.ft.
	MIR-H163-PA	MIR-H263-PA	MIR-154-PA	MIR-254-PA	MIR-554-PA	MLR-352H-PA
TEMPERATURE RANGE						
5°C above ambient to 80°C	■	■				
-10°C to 60°C			■	■	■	
5°C to 50°C and 10°C to 50°C						■
CONTROLLER						
Microprocessor	■	■	■	■	■	■
Setpoint Memory	■	■	■	■	■	■
PROGRAMMING						
Temperature	■	■	■	■	■	■
Humidity					■	■
Lighting					□	■
Program Memory	■	■	■	■	■	■
CONDITIONING						
Forced Air			■	■	■	■
Natural Convection (Fan Circulated)	■	■				
HUMIDIFICATION						
Ultrasonic Generated						■
LIGHTING						
Variable Intensity Fluorescent Lights						■
15 W Fluorescent Light			■	■	■	
DOOR						
Solid with Inner Glass Door	■	■				■
Solid with Viewing Window			■	■		
Solid with Viewing Window, with Cover					■	
CONSTRUCTION						
Stainless Steel Interior	■	■	■	■	■	■
Adjustable Sliding Shelves	■	■				
Adjustable Wire Shelves			■	■	■	■
Stackable			■			