

MVE Series




Features include


- Liquid sample storage
- Wide neck opening
- Lowest liftover height
- Largest LN₂ capacity
- Optional Battery Backup

Rack Layouts

4R9953 Layout



DF700 Layout



Featuring the TEC 3000

The MVE Series freezers provide stable cryogenic storage for up to 39,000 1.2 / 2.0 ml vials. These freezers provide maximum storage density and provide the industry's longest hold time. MVE Series are primarily designed for liquid storage.



www.chartbiomed.com
844-MVE-CRYO



	MVE 204 [†]	MVE 205	MVE 510	MVE 616	MVE 1426	MVE 1839												
Maximum Storage Capacity																		
1.2 & 2 ml Vials (Internally Threaded)	3,200	5,200	10,400	16,900	26,650	39,000												
Number of Large Racks 100 cell boxes	4	4	7	12	18	28												
Number of Mini Racks 25 cell boxes	-	-	4	4	10	8												
Number of Shelves per Rack	8	13	13	13	13	13												
High Security Straw Capacity (0.5 ml)	11,592	19,320	32,760	50,400	79,800	151,200												
Number of Canisters (73 mm)	23	23	39	60	95	150												
Goblets/Canister : Straws/Goblet	3:168	5:168	5:168	5:168	5:168	6:168												
Performance																		
LN ₂ Capacity w/o Inventory L	65	95	166	240	388	673												
Unit Dimensions																		
Neck Opening in. (mm)	16.00 (406)	16.00 (406)	20.70 (527)	25.10 (638)	31.75 (806)	39.40 (1002)												
Usable Internal Height in. (mm)	19.7 (502)	28.9 (735)	30.0 (762)	29.5 (749)	29.3 (745)	33.6 (854)												
Inner Diameter in. (mm)	16.00 (406)	16.00 (406)	20.70 (527)	25.1 (638)	31.75 (806)	39.40 (1002)												
Overall Height w/ AutoFill in. (mm)	30.8 (784) [†]	46.3 (1176)	46.7 (1185)	45.8 (1161)	46.8 (1188)	54.0 (1372)												
Door W x D Requirement** in. (mm)	18.0 (457)	23.6x22.7 (601x577)	30.6x29.7 (778x754)	33.5x34.6 (852x878)	40x41.6 (1017x1056)	46.8x48.9 (1187x1242)												
Weight Empty* lb. (kg)	79 (36)	195 (88)	281 (127)	320 (145)	490 (222)	750 (341)												
Weight Liquid Full* lb. (kg) est.	195 (88)	365 (166)	577 (262)	748 (339)	1181 (536)	1950 (885)												
Blood Bag Capacities																		
	Total Bags	Bags/ Frame	No. Frames	Total Bags	Bags/ Frame	No. Frames	Total Bags	Bags/ Frame	No. Frames	Total Bags	Bags/ Frame	No. Frames	Total Bags	Bags/ Frame	No. Frames	Total Bags	Bags/ Frame	No. Frames
791 OS/U Medsep (25 ml)	-	-	-	539	7	77	882	7	126	1,372	7	196	2,226	7	318	3,968	8	496
4R9951 (50 ml)	-	-	-	228	6	38	420	6	70	612	6	102	924	6	154	1,856	8	232
4R9953 (250 ml)	-	-	-	128	4	32	224	4	56	336	4	84	512	4	128	1,010	5	202
4R9955 (500 ml)	-	-	-	96	4	24	160	4	40	248	4	62	424	4	106	770	5	154
DF200 (200 ml)	-	-	-	80	4	20	128	4	32	200	4	50	336	4	84	610	5	122
DF700 (700 ml)	-	-	-	40	4	10	68	4	17	116	4	29	168	4	42	380	5	76


TWO Year Parts Warranty • FIVE Year Vacuum Warranty

Conforms to MDD 93/42/EEC, the Medical Device Directive for the EU. Freezer systems UL/C-UL Listed.

* Without inventory

**Minimum width required for vessel to pass through opening. Footprint may vary. Contact Tech Service for detailed drawings.

† Not available with Auto Fill.



Visit www.chartbiomed.com for more information.

© 2019 Chart Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.