

# Portable Anesthesia Machine



The leader in anesthesia equipment and supplies.

An anesthesia solution you can take anywhere.

Designed to be lightweight and easy-to-use, our Portable Anesthesia Machine (PAM) enables investigators to easily carry the device throughout the laboratory or into the field for special procedures. Our PAM is intended to be used with non-rebreathing systems, including the Universal Control Arm.

## Portable Anesthesia Machine

14.5" high x 11" wide x 8.5" deep  
07-8914722

Optional PAM accessories	Item #
Oxygen Flush Kit with 50 LPM Safety Restrictor	07-8910883
Dual Diverter Manifold	07-8914724
Tri Diverter Manifold	07-8914725
Nebulizer	07-8914735
Universal Control Arm Non-rebreathing Bain Block	07-8914012
Precision Vaporizers	07-8703584, 07-8703592, 07-8703600, 07-8703618
Multiplex Delivery Systems	07-8914732, 07-8914733, 07-8918185



PAM

***We can help you build your system.***

Call our experts today to place your order: 800.877.8989

