

Ultra Series

Ultra

The Ultra Basic Series is our value line of laboratory refrigerators and freezers for facilities that require quality equipment but face tighter budgets. Offering exceptional performance and temperature reliability, these units also feature:

- Microprocessor temperature controller, with digital temperature display
- Remote alarm contacts
- Forced air draft
- Remote alarm contacts
- Magnetic door gasket for positive seal
- 2" swivel casters
- LED interior lights
- Probe access port (3/8")
- UL/C-UL listed
- CFC-free refrigerant

Featured Display:



UltraElite

LHE

The Ultra Elite Series is a step above the Ultra Basic line. With a slightly different look and a larger temperature display, this series includes a majority of the same features as the Ultra Basic Series plus a few enhancements, such as:

- Includes Ultra Series features
- Sample simulated digital temperature display; .1° resolution
- Adjustable temperature display in celsius or fahrenheit
- Sensor error alarm
- Power failure alarm
- Resettable 1-point min/max history
- Alarm mute with ring back
- Swivel casters (glass door models only)
- HFC-free refrigerant (featured on glass door models only)

Featured Display:



UltraTouch

The Ultra Touch Series is comprised of our state-of-the-art laboratory refrigerators and freezers. These premium units provide the same performance and features as the Ultra Elite Series with the addition of:

- Includes Ultra & Ultra Elite Series features
- A touchscreen temperature display & alarm module
- Password-protected control panel
- Active temperature scale bar display
- Digital 24-hour temperature chart
- Temperature & alarm event data logging
- USB port for data transfer & software updates
- Real-time operational display clock
- 4-20mA output
- Adjustable alarm volume
- Door ajar alarm
- Low battery alarm

Featured Display:



Futura Series

Futura SILVER

The Futura Silver Series is the economic choice for all of one's laboratory needs. From undercounter refrigerators and freezers to hazardous material storage units, there is a model to suit everyone's needs both reliably and affordably. Diverse in features*, most models include:

- Digital Temperature Display
- High/low audible & visual alarms
- Remote alarm contacts
- Adjustable shelves
- Keyed door locks

*Features standard on most models



Futura LD

The Futura LD Series is comprised of full-featured laboratory refrigerators and auto defrost freezers that are ideal for the storage of temperature sensitive materials. All LD series models are ENERGY STAR certified and feature Hydrocarbon refrigerants.

- Microprocessor temperature controller
 - Battery backup included
- Alarm: Audible and visual high and low temperature alarms
- Remote alarm contacts
- Keyed Door Lock
- Energy Star high performance listed
- Hydrocarbon, natural refrigerant (R290)
- Forced draft air circulation
- Next generation monoblock compressor system (operating sound level 52 dB)
- LED interior lights
- Probe access port (1/2")



Futura PLUS+

The Futura Plus+ Series offers a high-end selection of refrigerators and freezers to suit a variety of laboratory needs. Accessories help to maximize storage capacity and efficiency while standard features* of this series include:

- Temperature controller with digital display
- High/low audible & visual alarms
- Remote alarm contacts
- Hydrocarbon, natural refrigerant (R290)
- Heavy-duty, swivel locking casters
- Adjustable shelves
- Interior lights
- Keyed door locks
- Door ajar alarm (hinged doors only)

*Features standard on most models

