


CENTRIFUGE END OF YEAR SALE

45% Off NuWind General Purpose Centrifuges and Accessories

Our overstock of centrifuges is your gain. Take advantage of this year end deal before time runs out.

Details: Offer Valid till December 31, 2019. Pricing valid until quoted expiration date.



MODELS AVAILABLE



HIGHLIGHTS

Max Speed	10,700 RPM
	13,000 xg
Max Capacity	40 x 50ml
	96 x 15ml
	156 x 10ml
	200 x 5ml

NU-C300V
4-Liter Benchtop Ventilated



HIGHLIGHTS

Max Speed	12,100 RPM
	16,660 xg
Max Capacity	40 x 50ml
	96 x 15ml
	156 x 10ml
	200 x 5ml


NU-C300R
4-Liter Benchtop Refrigerated



HIGHLIGHTS

Max Speed	12,100 RPM
	16,660 xg
Max Capacity	40 x 50ml
	96 x 15ml
	156 x 10ml
	200 x 5ml

NU-C300RF
4-Liter Floor Refrigerated



HIGHLIGHTS

Max Speed	15,300 RPM
	22,380 xg
Max Capacity	28 x 50ml
	56 x 15ml
	76 x 10ml
	92 x 5ml

NU-C200V
2-Liter Benchtop Ventilated



HIGHLIGHTS

Max Speed	18,000 RPM
	30,980 xg
Max Capacity	28 x 50ml
	56 x 15ml
	76 x 10ml
	92 x 5ml

NU-C200R
2-Liter Benchtop Refrigerated



POPULAR ACCESSORIES



Part #: NU-RX625
4 x 625 ml Swing Out Rotor

Accessories for NU-RX625



NU-B625
4 x 625 ml Swing Out Buckets (Qty 4)



NU-L625
4 x 625 ml Airtight Lids (Qty 4)



NU-BXDW
Microplate Carrier (Qty 2)

Inserts for NU-B625

Part Number	Description	Capacity
NU-2T50	5 x 50 ml Conical Bottom Tubes	20
NU-2T15C	11 x 15 ml Conical Bottom Tubes	44
NU-2T10	19 x 10 ml Round Bottom Tubes	76
NU-2T5	23 x 5/7 ml Round Bottom Tubes	92



Part #: NU-RA24-2
24 x 1.5 / 2 ml 45° Angular Rotor with Airtight Lid



Part #: NU-RA8-50
8 x 50 ml 25° Angular Rotor

REVIEWS

Kassa Beimnet
Hospital for Sick Children

Ease of use ★★★★★
After sales service ★★★★★
Value for money ★★★★★

Great product and admirable tech support!

Rating: 5.0 ★★★★★

Application Area: Specimen collection (bone marrow aspirate and peripheral blood) and DNA/RNA isolation procedures

"The product is very reliable/dependable and runs very smoothly. Despite having multiple features, it is very easy to use. I have already recommended it to others."

Review date: 06 Mar 2017 | Source: www.selectscience.net

Kim Ahrens
University of Florida

Ease of use ★★★★★
After sales service ★★★★★
Value for money ★★★★★

A great choice for any cell culture lab!

Rating: 5.0 ★★★★★

Application Area: Ficoll, general cell centrifuging flow cytometry.

"The NuWind bench top refrigerated centrifuge is a uniquely designed, quiet, user friendly addition to our lab. I love the capacity for 50 and 15 ml tubes, the touchscreen technology and the way it lets you know with a green check that your run has ended. It is programmable for different runs and the rotor easily changeable. I would definitely recommend this quality bench top centrifuge for the lab."

Review date: 13 Apr 2016 | Source: www.selectscience.net

Your Local NuAire Partner

LABRepCo

www.labrepco.com

